

CoolBox systems are bench top cooling workstations that provide sample cooling or freezing without ice, electricity or batteries! CoolBox systems are versatile and can accommodate a variety of sample formats and temperatures. The modular design enables the use of CoolRack® and CoolSink® thermo-conductive sample modules to hold microfuge tubes, cryogenic vials, PCR tubes and plates, assay plates and more. The internal cooling or freezing cartridges provide hours of 0.5 to 4°C cooling or sample freezing. CoolBox systems may also be used with dry ice to create a compact, portable, snap-freezing workstation for bacteria, virus and proteins.



Samples directly in ice

Non-reproducible method. Insulating air pockets create non-uniform temperatures. Unorganized and at risk of contamination.



CoolRack® module on ice

Reproducible method. Alloy rapidly adapts to ice temperature and holds all samples at the same temperature. No direct ice contact.



Ice-free CoolBox™ XT workstation

Reproducible method. Uniform and consistent sample temperature up to 16 hours with re-usable cooling core. No ice, no batteries, no electricity.

CoolBox Systems


| | CoolBox XT & 2XT | CoolBox MP | CoolBox 30 | XT Starter |
|------------------------------------|------------------|------------|------------|------------|
| Holds Tubes | • | • | • | • |
| Holds Plates | • | • | | • |
| o - 4°C cooling with lid closed | 16 hours | 10 hours | 10 hours | n/a |
| o - 4°C cooling with lid open | 10 hours | 5 hours | 4 hours | 4 hours |
| <0°C freezing with lid open | 5 hours | | 3 hours | |
| Details | Page 9 | Page 10 | Page 10 | Page 8 |

CoolBox Performance

Duration Test: Open Lid 20100 2 4 6 8 10 Time (hr)

CoolBox XT workstation performance with XT Cooling Core. Samples stay cold for over 10 hours with lid open and over 16 hours with lid closed.

How It Works



Each CoolBox system consists of a base, cooling cartridge/core and a lid. CoolBox XT workstations also include a magnetized collar and a 1 to 8°C temperature indicator for visual assurance of proper temperature. Each CoolBox is designed to receive thermo-conductive sample modules that sit directly on the cartridge/core keeping samples at desired temperature for hours. The cartridge/core is pre-frozen in a -20°C freezer.

Complete "Ready to Use" CoolBox Systems

Easily eliminate ice from your workflow by choosing one of these "ready to use" CoolBox ice-free cooling systems. Each system includes insulating box and lid, cooling cartridge or core, and sample module(s).



Ordering Information

Complete CoolBox XT and 2XT Workstations

| Application | For Use With | Description and Contents | Item No. Box 0 |
|--|--|--|----------------|
| Enzyme, PCR reagent or other biological specimen cooling | 24 x 1.5 mL or 2.0 mL microfuge tubes | MicroTube 24 CoolBox XT Workstation CoolBox XT System, SBS-compatible CoolRack XT M24 | BCS-576 |
| Cell thawing, tissue culture | 24 x 1.5 mL or 2.0 mL cryogenic vials (1.0 to 2.0 mL), with "gripping" wells for one-hand tube opening/closing | CryoTube 24 CoolBox XT Workstation CoolBox XT System, SBS-compatible CoolRack XT CFT24 | BCS-575 |
| PCR/qPCR plate preparation | One 96-well PCR plate, 12 strip wells or 96 0.2 mL tubes | PCR96 CoolBox XT Workstation CoolBox XT System, SBS-compatible CoolRack XT PCR96 | BCS-570 |
| PCR/qPCR strip well and reagent tube preparation | 6 strip wells and 12 microfuge tubes (1.5 to 2.0 mL) | PCR Strip CoolBox XT Workstation CoolBox XT System, SBS-compatible CoolRack XT M-PCR | BCS-572 |
| PCR/qPCR plate and reagent tube preparation | One 96-well PCR plate or strip wells or tubes and 24 x 1.5 mL or 2.0 mL microfuge tubes | PCR Sample Preparation CoolBox 2XT Workstation CoolBox 2XT System, SBS-compatible CoolRack XT M24 | BCS-573 |

Complete CoolBox MP Systems

| Application | For Use With | Description and Contents | Item No. Box Color |
|--|---|--|--------------------|
| PCR/qPCR plate preparation | One 96-well PCR plate or strip wells or tubes | PCR96 CoolBox System CoolBox MP System, SBS-compatible CoolRack XT PCR96 | BCS-145 |
| PCR/qPCR 384-well plate preparation | One 384-well PCR plate | PCR384 CoolBox System CoolBox MP System, SBS-compatible CoolRack XT PCR384 | BCS-146 |
| Tissue culture, cell staining, other assays | One 96-well or 384-well flat-bottom plate | 96F CoolBox System CoolBox MP System, SBS-compatible CoolSink XT 96F | BCS-147 |
| Bradford protein, other colorimetric or fluorescent assays | One 96-well u-bottom plate | 96U CoolBox System CoolBox MP System, SBS-compatible CoolSink XT 96U | BCS-148 |

Complete CoolBox 30 Systems

| Application | For Use With | Description and Contents | Item No. Box Color |
|--|--|---|--------------------|
| Enzyme, PCR reagent or other biological specimen cooling | 30 x 1.5 mL or 2.0 mL microfuge tubes | M3o CoolBox System CoolBox 3o System, CoolRack M3o | BCS-133 |
| Enzyme, PCR reagent or other biological specimen cooling (conical 1.5 mL tubes only) | 30 x 1.5 mL conical microfuge tubes | M30-PF CoolBox System CoolBox 30 System, CoolRack M30-PF | BCS-134 |
| Cell thawing, tissue culture | 30 x 1.0 mL or 2.0 mL cryogenic vials, with "gripping" wells for one-hand vial opening and closing | CFT30 CoolBox System CoolBox 30 System, CoolRack CFT30 | BCS-166 |

XT Starter



XT Starter is an open-platform cooler that accommodates most CoolRack and CoolSink modules. The low profile and small footprint make it ideal for use in the hood. Keeps samples cold (0.5 to 4.0°C) up to four hours. 1 to 8°C temperature indicator provides visual assurance of proper temperature. XT Starter includes one XT Holder and one XT Cooling Core. To extend the cooling duration, keep an additional XT Cooling Core in the freezer and replace the Core as needed.

XT Starter

Includes: (1) purple XT Holder and (1) XT Cooling Core. Accommodates one CoolRack or CoolSink sample module. (*Order separately, see page 11-14*)

 Item No.
 Description
 Dimensions (Lx W x H)

 BCS-504
 ● XT Starter, purple
 External: 16.3 x 12.1 x 6.1 cm

BCS-513 • XT Holder only, purple
BCS-511 XT Cooling Core

CoolBox XT and 2XT Workstations



CoolBox XT and CoolBox 2XT workstations keep tube or plate samples cold for over 16 hours with the lid on, and over 10 hours with the lid off. Use optional XT Freezing Core to maintain frozen samples for up to 8 hours. Dry ice may be used in place of the cores to create a compact snap freezing workstation. CoolBox XT Microfuge 24 Workstation (left).




CoolBox XT workstation
Holds one sample module



CoolBox 2XT workstation
Holds two sample modules

CoolBox XT Performance

| Temperature Range | Cooling Source | Open Lid | Closed Lid |
|-------------------|------------------|---------------|---------------|
| 0.5°C to 4°C | XT Cooling Core | Over 10 hours | Over 16 hours |
| -20°C to o°C | XT Freezing Core | Over 5 hours | Over 8 hours |
| -78°C | 200 mL Dry Ice | Over 4 hours | Over 5 hours |



"Design Your Own" Workstation

CoolBox XT Workstation

Includes one (1) XT Cooling Core for 0.5 to 4° C cooling. Holds one (1) sample module. (*Order separately from table on the right*)

| . , | , | , , |
|-----------|--|---------------------------------|
| Item No. | Description | Dimensions (L x W x H) |
| BCS-502 | CoolBox XT, purple | Internal: 14.0 x 9.4 x 9.9 cm* |
| BCS-502G | CoolBox XT, green | External: 20.0 x 16.0 x 15.0 cm |
| BCS-502O | CoolBox XT, orange | |
| BCS-502PK | CoolBox XT, pink | |

CoolBox 2XT Workstation

Includes two (2) XT Cooling Cores for 0.5 to 4°C cooling. Holds two (2) sample modules. (Order separately from table on the right)

| Item No. | Description | Dimensions (L x W x H) |
|-----------|---|----------------------------------|
| BCS-503 | CoolBox 2XT, purple | Internal: 21.0 x 14.3 x 13.4 cm* |
| BCS-503G | CoolBox 2XT, green | External: 26.5 x 20.0 x 15.0 cm |
| BCS-503O | CoolBox 2XT, orange | |
| BCS-503PK | CoolBox 2XT, pink | |

 $[*]Internal\ height\ of\ open\ space\ when\ core\ is\ in\ the\ base.$

Sample Modules for CoolBox XT and CoolBox 2XT Workstation

| Item No. | Description | For Use With | Dimensions (L x W x H) |
|----------|-------------------|---|------------------------|
| BCS-535 | CoolRack XT M24 | 24 x 1.5 mL microfuge tubes | 12.8 x 8.5 x 3.8 cm |
| BCS-534 | CoolRack XT CFT24 | 24 x 1.0 or 2.0 mL cryogenic vials | 12.8 x 8.5 x 3.8 cm |
| BCS-F20 | CoolRack XT PCRo6 | one of-well PCR plate strip wells tubes | 127 × 8 6 × 2 F cm |

Note: All XT modules are SBS-compatible and may be placed on robotic instruments for brief cooling outside of CoolBox.

CoolRack XT PCR96 one 96-well PCR plate, strip wells, tubes BCS-523 CoolRack XT M-PCR 12 x 1.5 or 2.0 mL microfuge tubes /6 strip wells 12.7 x 8.6 x 3.8 cm BCS-538 CoolRack XT PCR384 one 384-well plate 12.7 x 8.6 x 1.9 cm 12 injectable cell therapy ampules BCS-265 CoolRack SV10 12.7 x 8.6 x 3.8 cm one 96 or 384-well flat-bottom plate BCS-536 CoolSink XT 96F 12.8 x 8.5 x 1.7 cm BCS-537 CoolSink XT 96U one 96-well u-bottom plate 12.8 x 8.5 x 1.7 cm

Tall Tube Modules for CoolBox XT and CoolBox 2XT Workstation

Note: Tall modules require an XT Extension Collar for lid closing. See next page for details

| Twose. Tall Modules require all XT Extension Collar for his closing. See next page for details. | | | |
|---|---------------|--|----------------------|
| Item No. | Description | For Use With | Dimensions (Lx Wx H) |
| BCS-232 | CoolRack L* | 12 centrifuge tubes (15 mL), insulated | 13.7 x 9.5 x 11.8 cm |
| BCS-153 | CoolRack 15ml | 9 centrifuge tubes (15 mL) | 8.9 x 8.9 x 10.7 cm |
| BCS-154 | CoolRack 50ml | 4 centrifuge tubes (50 mL) | 8.9 x 8.9 x 10.7 cm |
| BCS-235 | CoolRack LV* | 12 blood tubes (13x75 mm), insulated | 13.7 x 9.5 x 9.6 cm |
| BCS-157 | CoolRack VS13 | 9 blood tubes (13x75 mm) | 8.9 x 8.9 x 6.1 cm |
| BCS-155 | CoolRack V13 | 9 blood tubes (13x100 mm) | 8.9 x 8.9 x 8.4 cm |
| BCS-156 | CoolRack V16 | 9 blood tubes (16x100 mm) | 8.9 x 8.9 x 8.4 cm |

^{*}No XT Extension Collar necessary unless lid closing is required.

CoolBox XT and CoolBox 2XT Optional Accessories



CoolBox XT Extension Collar

For use with CoolBox XT when using tall tube modules. Dimensions (L x W x H): Internal: 14.0 x 9.4 x 7.0 cm

External: 20 0 x 16 0 x 7 0 c

| External: 20.0 x 16.0 x 7.0 cm | | |
|--------------------------------|--|--|
| Item No. | Description | |
| BCS-502-C | Extension Collar, purple | |
| BCS-502-CG | Extension Collar, green | |
| BCS-502-CO | Extension Collar, orange | |
| BCS-502-CPK | Extension Collar, pink | |



CoolBox 2XT Extension Collar

For use with CoolBox 2XT when using tall tube modules. Dimensions (L x W x H): Internal: $21.0 \times 14.3 \times 7.0 \text{ cm}$

External: 26.5 x 20.0 x 7.0 cm

| External. 20.5 x 20.0 x 7.0 cm | | |
|--------------------------------|--|--|
| Item No. | Description | |
| BCS-503-C | Extension Collar, purple | |
| BCS-503-CG | Extension Collar, green | |
| BCS-503-CO | Extension Collar, orange | |
| BCS-503-CPK | Extension Collar, pink | |
| | | |



XT Cores


| Item No. | Description |
|----------|---------------------------------|
| BCS-511 | XT Cooling Core, For 0.5 to 4°C |
| BCS-512 | XT Freezing Core, For below o°C |

CoolBox MP Systems



CoolBox MP system keeps sample tubes or plates cold (0.5 to 4.0°C) for approximately 10 hours with the lid closed, and approximately 5 hours with the lid open. To extend the cooling duration, keep an additional cooling cartridge in the freezer and replace the cartridge as needed.





CoolBox MP Performance

| Temperature Range | Cooling Source | Open Lid | Closed Lid |
|-------------------|-------------------|---------------|----------------|
| o.5°C to 4°C | Cooling Cartridge | Up to 5 hours | Up to 10 hours |
| -78°C | 200 mL Dry Ice | Up to 4 hours | Up to 5 hours |

"Design Your Own" System

CoolBox MP System

Includes: CoolBox MP base and lid, one blue cooling cartridge. Accommodates one sample module (Order separately from table on the right.)

| | 13 | 0 / |
|-----------|--|-----------------------------------|
| Item No. | Description | Dimensions (L x W x H) |
| BCS-144 | CoolBox MP, purple | Internal: 13.3 x 8.8 x 7.3 cm* |
| BCS-144G | CoolBox MP, green | External: 20.3 x 15.25 x 15.25 cm |
| BCS-144O | CoolBox MP, orange | |
| BCS-144PK | CoolBox MP, pink | |
| | | |

^{*}Internal height of open space when core is in the base.

Cooling Cartridge

| item No | Desc | rip | lion | | | | |
|---------|------|-----|--------|-----|------|-------|------|
| | | | | | | | |
| BCS-152 | Cool | inσ | cartri | 100 | hlue | 1 ner | nack |
| | | | | | | | |



Sample Modules for CoolBox MP System

| Item No. | Description | For Use With | Dimensions (Lx W x H) | |
|--|---------------------|--|-----------------------|--|
| BCS-125 | CoolRack M15 | 15 x 1.5 or 2.0 mL microfuge tubes | 10.2 x 6.4 x 3.8 cm | |
| BCS-127 | CoolRack M15-PF | 15 conical 1.5 mL microfuge tubes | 10.2 x 6.4 x 3.8 cm | |
| BCS-535 | CoolRack XT M24* | 24 1.5 mL microfuge tubes | 12.8 x 8.5 x 3.8 cm | |
| BCS-126 | CoolRack CF15 | 15 cryogenic vials | 10.2 x 6.4 x 3.8 cm | |
| BCS-534 | CoolRack XT CFT24* | 24 cryogenic vials, with "gripping" wells | 12.8 x 8.5 x 3.8 cm | |
| BCS-529 | CoolRack XT PCR96* | one 96-well PCR plate or 12 strip wells | 12.7 x 8.6 x 2.5 cm | |
| BCS-523 | CoolRack XT M-PCR* | 12 x 1.5 or 2.0 mL microfuge tubes and 6 strip wells | 12.7 x 8.6 x 3.8 cm | |
| BCS-538 | CoolRack XT PCR384* | one 384-well PCR plate | 12.7 x 8.6 x 1.9 cm | |
| BCS-536 | CoolSink XT 96F* | one 96-well or 384-well flat-bottom plate | 12.8 x 8.5 x 1.7 cm | |
| BCS-537 | CoolSink XT 96U* | one 96-well u-bottom plate | 12.8 x 8.5 x 1.7 cm | |
| BCS-231 | CoolRack 96xo.5ml | 0.5 mL 2D storage tubes | 13.2 x 8.9 x 3.6 cm | |
| BCS-149 | CoolRack 96x1ml | 1.4 mL 2D storage tubes | 13.2 x 8.9 x 4.3 cm | |
| *All XT modules are SRS compatible and may be placed an robotic instruments for brief coaling outside of CoalRay | | | | |

^{*}All XT modules are SBS-compatible and may be placed on robotic instruments for brief cooling outside of CoolBox.

CoolBox 30 Systems




CoolBox 30 System keeps sample tubes cold (0.5 to 4.0°C) for approximately 4 to 10 hours. Use the optional Freezing Cartridge to maintain frozen samples below 0°C for approximately 3 to 6 hours.

CoolBox 30 Performance

| Temperature Range | Cooling Source | Open Lid | Closed Lid |
|-------------------|--------------------|---------------|----------------|
| o.5°C to 4°C | Cooling Cartridge | Up to 4 hours | Up to 10 hours |
| -20°C to 0°C | Freezing Cartridge | Up to 3 hours | Up to 6 hours |
| -78°C | 200 mL Dry Ice | Up to 4 hours | Up to 5 hours |





"Design Your Own" System

CoolBox 30 System

Includes: CoolBox 30 base and lid, blue cooling cartridge. Accommodates 15-well or 30-well CoolRack tube modules. (Order separately from table on the right.)

| , | 1 3 | 0 / |
|-----------|--|----------------------------------|
| Item No. | Description | Dimensions (L x W x H) |
| BCS-130 | CoolBox 30, purple | Internal: 12 x 10.4 x 7.9 cm* |
| BCS-130G | CoolBox 30, green | External: 16.5 x 15.25 x 12.7 cm |
| BCS-130O | CoolBox 30, orange | |
| BCS-130PK | CoolBox 30, pink | |

 $[*]Internal\ height\ of\ open\ space\ when\ core\ is\ in\ the\ base.$

Cooling Cartridge

For 0.5 to 4°C cooling
Item No. Description
BCS-132 Cooling cartridge, blue, 3 per pack



Sample Modules for CoolBox 30 System

| Sample Modules for Coolbox 30 System | | | | |
|--------------------------------------|--------------------------|--|----------------------|--|
| Item No. | Description | For Use With | Dimensions (LxWxH) | |
| BCS-125 | CoolRack M15 | 15 x 1.5 mL or 2.0 mL microfuge tubes | 10.2 x 6.4 x 3.8 cm | |
| BCS-127 | CoolRack M15-PF | 15 x conical 1.5 mL microfuge tubes | 10.2 x 6.4 x 3.8 cm | |
| BCS-108 | CoolRack M30 | 30 x 1.5 mL or 2.0 mL microfuge tubes | 12.0 x 10.2 x 3.8 cm | |
| BCS-128 | CoolRack M30-PF | 30 x conical 1.5 mL microfuge tubes | 12.0 x 10.2 x 3.8 cm | |
| BCS-137 | CoolRack M30-PF (500 uL) | 30 x conical 500 uL tubes | 12.0 x 10.2 x 3.8 cm | |
| BCS-126 | CoolRack CF15 | 15 x 1.5 mL or 2.0 mL cryogenic vials or 12x75 FACS tubes | 10.2 x 6.4 x 3.8 cm | |
| BCS-138 | CoolRack CFT30 | 30 x 1.5 mL or 2.0 mL cryogenic vials or 12x75 FACS tubes, with "gripping" wells | 12.0 x 10.2 x 3.8 cm | |

Freezing Cartridge

For maintaining frozen samples below o°C

| Item No. | Description |
|----------|---------------------------------------|
| BCS-131 | Freezing cartridge, green, 3 per pack |

