


Sliding Door Pharmaceutical Refrigerators

MPR Pharmaceutical Refrigerators

MPR-S150H-PE | MPR-S300H-PE



#### Introduction of PHCbi New MPR Series

Adoption of HC refrigerants and inverter compressor:

New lineup of MPR based on high energy-saving performance

(compared to existing models)









MPR Pharmaceutical Refrigerator with Freezer

#### Introduction of MPR-S150H/S300H

# Adoption of HC refrigerants and inverter compressor High energy-saving performance (compared to the existing models) Reduce condensation on the glass

1

HC refrigerant consumes about 64% less energy than conventional models \*Compared to MPR-S313-PE AT:23 °C SV: 5 °C 230V / 50 Hz

By HC refrigerants and inverter compressor 64% energy-saving (compared with existing models). Management of Freon gas is unnecessary.

2

#### Reduce condensation on the glass

Use of highly insulated argon gas filled pair glass and air flow from the bottom of the door significantly reduces condensation on the door glass. Ensure high temp. performance.



#### Improved usability

New display with organic EL has good visibility and intuitive operation. USB port and log functions simplified temperature data management.





PR Pharmaceutical

MPR-S150H-PE [165L] MPR-S300H-PE [345L]



### **Environmental compliance and energy saving**

Adoption of HC refrigerants and inverter compressor (Refrigerator compartment)

# Approx. 64% reduction compared to the existing model

\*1 : Compared to MPR-S313-PE \*2 : AT:23 °C SV: 5 °C 230V / 50 Hz



HC refrigerants minimize environmental impact without cooling performance degradation



The newly equipped inverter compressor can control the speed of the compressor to reduce power consumption by reduction

the on/off cycle which consumes power.




















#### Improvement of Temp. control performance

**Inverter control + New platform adoption** 

## Ensure 2 to 8°C with Peak to Peak

Reduce Temp. fluctuation by adopting inverter compressor Improved Temp. distribution by adopting new platform


Cabinet
Peak to Peak



#### MPR-S300H-PE


Peak to Peak: 4.4°C Max 7.2°C, Min 2.8°C

Adjusts the compressor rotation speed in 5 steps and suppresses the peak temperature



#### MPR-S313-PE

Peak to Peak: 5.9°C Max 8.3°C, Min 2.4°C





(Condition: AT23 °C, SV5 °C, 9 point method)

рнсы


Copyright © 2019 PHC Corporation All Rights Reserved.

## Temperature stability for medication & sample storage

Defrosting operation is possible without turning off the power even during defrosting, stabilizing the inner temperature.

Cycle defrost


Power consumption is further reduced by detecting the temperature of the evaporator and providing time to turn off both heater and compressor.





## Reduced condensation on glass door

Use of highly insulated argon gas filled two layer glass and air flow from the bottom of the door

Significantly reduces condensation on the glass







### **Improved Usability**

# New display with organic EL

More detailed temperature display with increment of 0.1°C



#### **Improved Usability**

# New display with organic EL

More detailed temperature display with increment of 0.1°C



- Flexible alarm settings by 1°C in high temp. / low temp.
- Display of alarm message
- Display of simplified temp. graph







#### **Improved Usability**

# **New display with organic EL**

More detailed temperature display with increment of 0.1°C



**USB** port and data log function

- Function that logs temperature, alarm, door opening/closing
- Can store a max. 3 months of data (1min. interval)
- Can export data as a CSV file



Temperature logs (CSV)

Display of

**Alarm logs** 

08:45~ W10 ▼

Display of Door Opening/closing logs

Alarm 18/03/05 12:00~ A01 Door Open 18/01/25 05:30~ W03 Ref.SNSR Open 18/01/01 08:45~ W10 ▼ Power Failure

## New platform improved storage efficiency


Increased capacity while maintaining an installation area and a height that is easy to access






## **New platform improved storage efficiency**





## PHCbi MPR series line-up







**Life Science Innovator Since 1966**