

# **CLIMACELL** EVO

Climatic chamber with forced air circulation, cooling and controlled humidity









# Tradition, quality,

#### innovation

Since its establishment in 1921, BMT Medical Technology s.r.o., the traditional manufacturer of medical and laboratory technology, has been gradually transformed from a small regional company to an international corporation.

In 1992, it became a member of the European MMM Group which has been operating on the world markets since 1954 as an important supplier of systems for the health care industry, science and research. With its comprehensive offer of products and services, sterilization and disinfection devices for hospitals, scientific institutes, laboratories and pharmaceutical industry, MMM Group has established itself as an outstanding quality and innovations producer on the global markets.

The knowledge and experience gained during the implementations of individual supplies for our customers all over the world, and the technical innovations have been permanently and positively influencing the development, construction and production of our devices. High level of our work has also been confirmed by the number of patents and industrial designs as well as an easy implementation of individual device adjustments.

MMM Group excellence in medical and laboratory technology.

# **CLIMACELL** EVO

### Climatic chamber with a wide range of applications

Air-conditioned chambers of the CLIMACELL line provide all and any conditions for exact and reproducible simulation of various climatic conditions. Thanks to a wide range of adjustable parameters -20°C up to +100°C- of temperature and 10-95% of humidity, possibility of CO2 regulation and sterilisation at 160°C and numerous variety of options and accessories, the new generation of the climatic chamber Climacell Evo becomes an ideal tool for simulation of conditions in many fields of human activities. Simple control via touch screen, precise regulation and wide possibilities of data outputs meet the most demanding conditions of pharmaceutical industry and they also allow user friendly simulation of simple requirements for plants growing. The devices offer an interesting alternative for expensive testing chambers and testing rooms. The microprocessor-controlled system of humidification and dehumidification together with high-performance programmable system of exposition lighting quarantee excellent homogenous parameters for tests and growth conditions.

Meeting the requirements of regulations 2006/95/EC, 2004/108/EC, ICH 279/95 Option 2, FDA 21 part 11.



### Applications



Pharmaceutical Industry

Stability testing and photo stability testing according to ICH 279/95 option 2, long term storage



**Electronic Industry** 

Testing of durability of electronic boards and printed circuits



**Cosmetic Industry** 

Durability testing, testing of cosmetic products or primary materials stability



**Automotive Industry** 

Testing of materials ageing - tyres, sealing, etc.



**Construction Industry** 

Long-term testing of quality and ageing of materials in construction industry - cement, paints, asphalt, construction plastics, glues, etc.



Zoology

Simulation of conditions for sea organisms research - seaweed or cultivation of insect eggs.



General and Applied Industry (research field)

E.g. cultivation of tissue cultures human or animal ones



Studies of germination, green plants growing for further research



Food and Beverage Industry

Testing of food quality under simulated transport or storage conditions - export of fruits, etc.




Field of Metrology and Quality in Industry Checking and calibration of

industrial measuring gauges

**Packing Material Industry** 

Long-term testing of packing technologies






### The new control system offers

- 5,7 inch (14,5 cm) touch screen display
- Microprocessor fuzzy logic process control
- Intuitive control via colour icons
- Graphic configuration of a new program
- Transparent displaying of data course at the cycle
- Protective thermostat class 3
- · Acoustic and visual alarm
- Multi-level users administration (corresponding to FDA 21 Part 11)
- · Keyboard lock against unauthorised handling

- Data encryption and non-manipulability (corresponding to FDA 21 Part 11)
- Up to 100 programs and up to 100 segments for each program
- 30 day data logger in graphic and numeric form
- On-line or off line data export
- Prepared service programs for fast diagnostics of faults
- Easy service diagnostics including remote
- Multi-language communication
- Direct printing of protocols in pdf format

- Easy user configuration of the device
- SD memory card, USB Host and RS 232 standardly included
- WIFI connection, USB device or Ethernet interface with own IP address for remote data transfer, control and diagnostics (optional equipment)
- · Programming of ramps, real time and cycling
- Fan setting 0-100%
- Main ON/OFF switch for security reasons
- Device functionality LED indicator

### Connectivity



**Connection of software** Warmcomm 4.0 **Smartphones** remote (internet)



# **Data Outputs**

Thanks to the most up-to-date components of electronic, the device Climacell EVO does not have any data peripherals connection limit. The basic configuration contains traditional and reliable RS 232, USB Host and the SD card as the data carrier. The device can be easily extended with the Wifi 802.11b/g module with up to 100 meters reach, there is also available the USB device for bi-directional USB communication and for remote connection there is the Ethernet (RJ 45) connection. Proper IP address allows easy connection to PC or selected WIFI printer, respectively other usual data periphery (Smartphone, Netbook, etc.). Thanks to the open platform and adjusted data format it is also possible to configure remote connection and to work with on-line data in remote mode (internet).

# WARMCOMM 4.0

printer via RS 232

Integrated data Compatible thermo-

Universal Data Administration with devices of the MMM Group

Specified desk-top

printer (USB/WiFi)



logger, SD card











- connectible to all the devices of the MMM Group
- stable platform of the SQL library
- user-friendly environment
- connection of up to 25 devices
- bilateral communication data monitoring and device control
- compatibility with older lines of heat technology devices
- Client-Server architecture
- service module for local and remote diagnostics
- three levels depending on client's requirements (B-P-F)
- compliance with FDA CFR 21 Part 11 (version F)
- web support, on-line updates
- protected licence policy
- unpretentious to HW requirements, compatible with MS Windows and UNIX



# **CLIMACELL** EVO

#### Comfort machine with


### superior parameters

MMM Group offers traditionally fully ranged size of the cabinet, from personal size 111 litres, up to new size 1212 litres, with the best ratio cost/performance. Patented vertical air flow with preheating chamber and asymmetrically perforated panels ensure the well proven vertical spiralled air flow with the best spatial homogeneity.

Deep experience of the factory engineers and many years of careful development help with sophisticated fuzzy logic control system. By means of the fuzzy logic are continually evaluated the current process conditions as size of chamber, set parameters, quantity of the samples inside and herewith optimizing heating, cooling and steaming performance.

High pressure steam generator in new easy accessible position and newly designed powerful freezing coil regulate the relative humidity quickly in full range from 10–95 % RH, according the customer set, and without significant temperature interference.

Practical large and popular door handle, robust wheels with brakes and 170° openable main door(s) contributes to high user friendly character of the device. Light grey with light blue device colours highlighted by dark blue smiley control panel cause a pleasant feeling of harmony in the user every morning













0

# **Programmable lighting**

The new generation of climatic chamber Climacell EVO offers wide possibilities of light options applications. Together with control of temperature, humidity, ventilation or CO2 it is possible to control and program even intensity of visible or ultra-violet exposure in the form of LED diodes.

#### Light shelves

In case of testing the photo-stability for pharmaceutical industry, the most suitable option is the use of light shelves for exposition of samples. The lighting is in compliance with standards ICH Q1A and Q1B, option 2 and relevant European standards. There are available light shelves with visible light, ultra-violet light and their combinations. It is possible to set different temperature or humidity for each segment and accordingly, it is possible to set different level of lighting for each segment. That allows programming of cycles of day and night simulations with gradual light intensity increase and decrease.







### **Exposition Lighting in Doors**

For applications with necessary equal lighting of the chamber, for example for plants growth, there is available an option with lighting in the device door. Just like in the light shelves, even exposition light in door allows separate regulation for each segment within the range of 0-100% in steps of 10%.




### **Light Sensors**


To meet the most strict requirements for light exposition measurability there is available the option of light sensors. One type of sensors is able to exactly measure the quantity of visible and the other the quantity of ultra-violet light. In case of such sensors installation it is possible to follow exact exposition doses received by samples on a display, in graphs or in the software Warmcomm 4.




### Light Spectrums Available

As there extend the possibilities of LED diodes manufacturers, there also extend the possibilities of their use in Climacell EVO. Together with usual possibilities of visible (white) light or ultra violet light we will try our best to get even LED diodes meeting your requirements towards other light colours. Please do not hesitate to contact us in such cases and we will be glad to discuss the possibilities with you.







# **Humidity Control**

**Steam Generator** 

of relative humidity.

The device allows steam generation and

its precise dosing to the chamber. Thanks

to our long-term experience in the field

of steam sterilisation we succeeded to

develop pressure steam generator able to

increase relative humidity in the chamber

in a precise, reliable and fast way. Steam

reservoir using the heating element. Then,

the chamber. The technology eliminates the

overshootings while reaching required level

overpressure is generated in the water

Climacell EVO is a climatic chamber – i.e. a device that is able to exactly and quickly regulate the quantity of humidity in the chamber. This is possible thanks to strong system of active increase and decrease of humidity in connection with the system of water supply to the device.

Unlike many other manufacturers we are using the separate cooling snake of the cooling system. The control system of humidity on the cooling coil condensates the valve releases exact volume of steam to



### **Humidity Reduction**

not engaged only in humidity increase, but we also focus on active humidity decrease, Climacell EVO is able to reduce humidity in the chamber using the cooling system while keeping a nearly constant temperature. The



#### Water Intake and Use

and condensed water is taken

where the pump releases it to

the water waste. Thanks to the

quickly.

efficient system Climacell EVO reaches

even the low relative humidity values very


back to the steam generator

In order to create the humidity exactly and reliably in the long term, the steam generator of Climacell EVO operates only with demineralized water. The access to such water can be solved in two ways. A standard solution means that you pour demineralized water to a barrel, delivered with each Climacell EVO and you connect the pump from the barrel to the connector on the rear side of the device. The other possibility includes connection of demineralized water intake from the laboratory water distribution system to the steam generator of Climacell EVO via reduction pressure valve. In both cases, the device automatically takes exact quantity of water as needed for humidity creation in the steam generator.



### Restrictions of temperature and relative humidity setting combinations

10 20 30 40 50 60 70 80 90 100 T (°C)



#### Accessories included

Each Climacell EVO is supplied with standard equipment which does not have to be additionally ordered and it makes a standard part of delivery:



Touch screen



Communication ports RS 232 and







Water barrel for distilled water



Reliable RH sensor



Multi-conductor temperature sensors



Sealing inner glass door



2 stainless steel trays

## **Optional Equipment**

Thanks to modular construction of our devices even Climacell EVO may be additionally equipped according to your preferences with many additional options. Climacell EVO may then serve as a chamber for testing of photo-stability, light simulation of day and night, processes with CO2 control, hot-air sterilisation, etc.

- 1. Hot-air sterilisation 160 °C
- 2. Additional cooling -20°C
- 3. Flexible temperature sensors
- 4. LED light shelves
- 5. Exposition lighting in doors
- 6. Light sensors of exposition

- 7. Defrosting system
- 8. CO2 control
- 9. Software Warmcomm 4
- 10. Data module USB device, Ethernet, wi-fi
- 11. Mechanic door lock
- 12. Electromagnetic door lock
- 13. Trays or shelves
  - 14. Access port 25/50/100 mm
  - 15. Programmable inner socket
  - 16. External printer
  - 17. Multi-point temp. / humidity measuring
  - 18. IQ/OQ protocols








































# **Technical parameters**



	CLIMACELL E\	/0 111, 20	2, 404, 70	7, 1212			
Technical data Internal space - chamber, stainless steel DIN 1.4301 (AISI 304)	Volume	approx. litres	111	222	404	707	1212
	Width	mm	540	540	540	940	3x540
	Depth	mm	370	520	520	520	520
	Height	mm	530	760	1410	1410	1410
Steam space volume		approx. litres	167	305	530	878	1781
Trays of stainless steel *)	Racks	max. No.	7	10	19	19	3x19
	Standard equipment	pcs	2	2	2	2	6
	min. distance between trays	Cm	7	7	7	7	7
Maximal allowed load of trays *)	Tray loading capacity	max. kg/ tray	20	30	30	50	30
	Loading capacity in total	kg/case	50	70	100	130	300
Number of doors	External solid / internal glass	pcs	1/1	1/1	1/1	2/2	3/3
External dimensions (including door, handle and wheels)	Width	mm	780	780	1100	1500	2630
	Depth	mm	755	885	885	885	885
	Height	mm	1177	1452	1888	1888	1910
Weight CLC 0°C	Net	approx. kg	110	140	240	280	-
	Gross	approx. kg	140	177	280	326	-
Weight CLC -20°C	Net	approx. kg	120	150	250	290	340
	Gross	approx. kg	150	187	290	336	395
Electric parameters	max. input CLC 0 °C	W	2050	2100	3150	3400	-
	max. input CLC -20 °C	W	1630	1780	2115	2640	3215
	mains 50/60 Hz	V	230	230	230	230	230
	Protective system		IP20	IP20	IP20	IP20	IP20
Temperature data Operation temperature	from 0.0 °C from -20.0 °C	to °C to °C	100 (decontamination 160°C) 100 (decontamination 160°C)				- -20 - +70
Temperature accuracy	In - Space at 10 °C	approx. (±) °C	<0,5	<0,5	<1	<1	<1,5
	at 37 °C	approx. ( <u>+</u> ) °C	<0,5	<0,5	<1	<1	<1,5
	In - Time	approx. ( <u>+</u> ) °C	<0,2	<0,2	<0,3	<0,4	<0,8
Heating up time to 37 °C from the ambient temperature		min	20	25	26	27	30
Cooling/down time from 22 °C to 10 °C		min	<21	<21	<21	<21	<21
Recovery time after 30 s of door opening according to DIN 12 880	at 37 °C	min	4	4	4	4	4
	at 50 °C	min	4	4	4	4	4
Relative humidity CLC	range	0/0	10-95	10-95	10-95	10-95	10-95
Heat emission	at 37 °C	approx. W	70	97	123	148	200
Complete device noise level		dB	46	50	56	58	60
CO2 concentration		0/0	0,1 - 20			-	-
Required pressure CO2		bar/psi	0,3-0,7/5-10			-	-

Construction changes reserved

Note: All technical data are related to 22 °C ambient temperature and ± 10% voltage swing (if not specified). For other parameters see section Electric connections. There occurs temperature and humidity variation in case of regular turbidity removal during the operation.

\*) Approx. 50 % of the tray area can be filled in a way a uniform air circulation is enabled inside the chamber.

## Make acquaintance with our further offers ...

# Laboratory ovens and incubators



#### **FCOCFII**

The highly cost-effective heating oven series for simple drying processes

#### DUROCEL

Special- purpose drying ovens DUROCELL with highly resistant EPOLON coating

#### VACUCELL

Vacuum drying ovens

#### **STERICELL**

Intended for hot air sterilization of materials under specifiedtemperature and duration.

#### VENTICE

Drying ovens with forced air circulation

#### INCUCELL / INCUCELL V

Suitable for safe treatment of microbiological cultures

#### **FRIOCELL**

Cooling incubators

#### CLIMACELL

Climatic chambers

#### CO2CELL

CO<sub>2</sub> atmosphere

## Sterilization and depyrogenation




#### VENTICELL IL

series of modular large-sized laboratory devices with the chamber volume of from 700 to 2000 liters. The devices are used for items sterilization at the temperature of up to 180  $^{\circ}$ C, or for items depyrogenation at the temperature of up to 300  $^{\circ}$ C and optional time mode. The devices can be used in laboratories, industry, pharmacy, and research.

## Steam sterilizers (autoclaves)



STERILAB
Small steam sterilizer, 25 l



UNISTERI HP IL
Medium-sized steam sterilizers, 73–254 l



STERIVAP HP IL Large steam sterilizers, 148–1490 I

