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Originalversion

OPERATING MANUAL

MODE D'EMPLOI

MANUAL DE INSTRUCCIONES

TITRONIC® 300

KOLBENBÜRETTE
PISTON BURETTE
BURETTE A PISTON
BURETA DE EMBOLO

SI Analytics

a xylem brand

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Operating Manual Page 53 .. 104

Important notes:

The operating manual is part of the product. Before initial operation, please carefully read and observe the operating manual and keep it. For safety reasons the product may only be used for the purposes described in these present operating manual. Please also consider the operating manuals for the devices to be connected

All specifications in this operating manual are guidance values which are valid at the time of printing. However, for technical or commercial reasons or in the necessity to comply with the statutory stipulations of various countries, the manufacturer may perform additions to the product without changing the described properties.

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TABLE OF CONTENT

1	Technical Specifications of the TITRONIC® 300 Piston Burette	55
1.1	Summary	55
1.2	Specifications Piston burette TITRONIC® 300	56
1.3	Warning and safety information	58
2	Unpacking and First Operation	59
2.1	Unpacking and First Operation of the Piston Burette.....	59
2.2	Connection and installing of Piston Burette and magnetic stirrer TM 50	59
2.3	Connecting the piston burette – Combination with Accessories and Additional Devices	60
2.3.1	Back panel of the TITRONIC® 300 piston burette	60
2.3.2	Connection ports of the TITRONIC® 300 piston burette	60
2.3.3	Connecting a printer	60
2.3.4	Connecting a USB device (manual controller, keyboard, memory device, hub)	60
2.3.5	Connection of analytical balances.....	60
2.4	Setting the Language of the Country	61
2.5	Dosing unit and Accessories	63
2.6	Installing the burette tip	64
2.7	Initial Filling or Rinsing of the Dosing unit	65
3	Working with the Piston burette TITRONIC® 300	67
3.1	Front Keyboard.....	67
3.2	Display	67
3.3	Manual controller.....	68
3.4	External PC Keyboard (Optional)	68
3.5	Menu Structure	69
3.6	Main Menu.....	71
3.6.1	Manual Titration.....	71
3.6.2	Dosage	73
4	Method Parameters	76
4.1	Method editing and new method	76
4.2	Default Methods	77
4.3	Copy Methods	77
4.4	Delete Methods	78
4.5	Change Method Parameters	78
4.5.1	Method type.....	78
4.5.2	Result (only for manual Titration)	79
4.5.3	Dosing parameters	84
4.5.4	Sample identification	85
4.5.5	Documentation	85
5	System Settings.....	87
5.1	Reagents/Dosing Unit	87
5.1.1	Dosing unit exchange.....	87
5.2	Global Memories	91
5.3	RS232 Settings	91
5.4	Connection of Printers.....	93
5.5	Date and Time.....	93
5.6	RESET	94
5.7	Device Information	94
5.8	System Sounds	95
5.9	Data exchange	95
5.10	Software Update.....	96

6	Connection of Analytical Balances.....	99
6.1	Balance data editor	99
7	Data Communication via RS-232- and USB-B interface.....	101
7.1	General Information.....	101
7.2	Chaining multiple devices — „Daisy Chain Concept“	101
7.3	Instruction Set for RS-Communication.....	102
8	Maintenance and Care of the TITRONIC® 300 Piston Burette.....	103
9	Storage and transportation.....	104
10	Recycling and Disposal	104

Notes to the Manual

The provided manual will allow you the proper and safe handling of the product.
For maximum security, observe the safety and warning instructions in the Instructions.

The pictogram  has the following meaning:

- Warning of a general danger to personnel and equipment.
- Non-compliance results (can result) in injury or material damage.

 Important information for device use.

 Refers to another part of the operating manual.

Status at time of printing

Advanced technology and the high quality of our products are guaranteed by a continuous development. This may result in differences between this operating manual and your product. We cannot exclude mistakes. We are sure you understand that no legal claims can be derived from the information, illustrations and descriptions.

A potentially more recent version of this manual is available on our internet website at www.si-analytics.com. The German version is the original version and binding in all specifications.

Copyright

© 2016, Xylem Analytics Germany GmbH.

Reprinting - even as excerpts - is only allowed with the explicit written authorization.

1 Technical Specifications of the TITRONIC® 300 Piston Burette

1.1 Summary

The TITRONIC® 300 is a piston burette and suitable for the following applications:

- Manual titration with or without calculation of the result
- Dosage
- It can be used as a dosing burette with the TitroLine® 7750 Titrator
- It can be used as a dosing or titrating burette in combination with the TitrSoft control software from version 3.1 on-going.

A variety of dosing and filling speeds can be set for each method.

Up to 3 user methods can be memorised in the device.

Solutions to be used:

Virtually, any liquids and solutions with a viscosity of $\leq 10 \text{ mm}^2 / \text{s}$ such as concentrated sulphuric acid may be used. However, one has to avoid the use of chemicals that may attack glass, PTFE or FEP or that are explosive, such as hydrofluoric acid, sodium azide or bromine! Suspensions containing high solids percentages may clog or even damage the dosing system.

General provisions:

The safety guidelines that are applicable to the handling of chemicals have to be observed under all circumstances. This applies in particular to inflammable and/or etching liquids.

Guarantee

We provide guarantee for the device described for two years from the date of purchase. This guarantee covers manufacturing faults being discovered within the mentioned period of two years. Claim under guarantee covers only the restoration of functionality, not any further claim for damages or financial loss. Improper handling/use or illegitimate opening of the device results in loss of the guarantee rights. The guarantee does not cover wear parts, as lobes, cylinders, valves and pipes including the thread connections and the titration tips. The breach of glass parts is also excluded. To ascertain the guarantee liability, please return the instrument and proof of purchase together with the date of purchase freight paid or prepaid.

1.2 Specifications Piston burette TITRONIC® 300

(Release: 30. August 2018)

CE Zeichen: As per the low voltage directive 2014/35/EU;
 test basis EN 61 010-1: 2011-07 for laboratory devices
 As per the EMC directive 2014/30/EU; test basis EN 61 326 Part 1: 2012
 As per the EMC directive 2011/65/EU; test basis EN 50 581: 2013-02
 FCC part 15B and ICES 003

Country of origin: Germany, Made in Germany

The following dissolvents/titration reagents are allowed to be used:

- All common titration solutions.
- As reagent water and all non-aggressive non-organic and organic fluids are allowed. If using combustible fluids fire please adhere to the Guidelines for Explosion Protection and Prevention of the chemical industry.
- For fluids with higher viscosity ($\geq 5 \text{ mm}^2/\text{s}$), lower boiling point or affinity to outgas, the filling and dosage speed can be adjusted.
- Fluids with viscosity over $20 \text{ mm}^2/\text{s}$ cannot be dosed.

Display: 3.5 inches -1/4 VGA TFT display with 320x240 pixels.

Power supply: power supply 100-240 V; 50-60 Hz,
 power input: 30 VA
 Input: 12 Volt DC, 500 mA.

Use the Power supply TZ 1853, Type No.: FW 7362M/12 only!

RS-232-C Interface: RS-232-C interface separated, Daisy Chain function available.
 Data bits: adjustable, 7 or **8** Bit (default: 8 Bit)
 Stop bit: adjustable, **1** or 2 Bit (default: 1 Bit)
 Start bit: static **1** Bit
 Parity: adjustable: even / odd / **none**
 Baud rate: adjustable: 1200, 2400, **4800**, 9600, 19200 (Default 4800 baud)
 Address: adjustable, (0 to 15, default: 01)

RS-232-1 for computer, input Daisy Chain
 RS-232-2 devices of SI Analytics, titrator TitroLine® 7750,
 - Burettes TITRONIC® 500, TITRONIC® 110 *plus*, TITRONIC® *universal/300*,
 - Balances of the types Mettler, Sartorius, Kern, Ohaus
 (for more, please contact SI Analytics)
 - Exit Daisy Chain

USB Interface: 1 x USB-type A and 1 x USB-type B

USB –Typ B for connecting a PC

USB –Typ A („Host“)
 for connecting:
 - USB keyboard
 - USB printer
 - USB manual controller,
 - USB data media e.g. USB stick
 - USB Hub to connect more USB devices

Stirrer: Plug connection with integrated low-voltage supply (15 V \approx) for the TM 50 Stirrer incorporated in the bottom of the casing of the TITRONIC® 300 piston burette

Housing material: Polypropylene
Front keyboard: Polyester
Housing dimensions: 135 x 310 x 205 mm (w x h x d), height incl. cylinder, excluding stirrer
Weight: ca. 2 kg for basic unit

Ambient conditions: Ambient temperature: + 10 ... + 40 °C for operation and storage
Humidity according to EN 61 010, Part 1:
Max. relative humidity 80 % for temperatures up to 31 °C,
linear decrease down to 50 % relative humidity at a temperature of 40 °C

Dosing units

Cylinder: Piston burette TITRONIC® 300 : 20 ml DURAN® (borosilicate glass 3.3) cylinder
Piston burette TITRONIC® 300 : 50 ml DURAN® (borosilicate glass 3.3) cylinder
UV protection: cover made out of TROGAAMID, blue, transparent

Valve: volume neutral cone valve made from fluorocarbon polymers (PTFE), TZ 3000

Hoses: FEP hose set, blue

Dosing accuracy: after DIN EN ISO 8655, part 3
Accuracy: 0.15 %
Precision: 0.05 %

1.3 Warning and safety information

The Piston Burette TITRONIC® 300 corresponds to protection class III. It was manufactured and tested according to DIN EN 61 010, Part 1, Protective Measures for electronic measurement devices and control devices and has left the factory in an impeccable condition as concerns safety technology. In order to maintain this condition and to ensure safe operation, the user should observe the notes and warning information contained in the present operating instructions. Development and production is done within a system which meets the requirements laid down in the DIN EN ISO 9001 standard.

For reasons of safety, the Piston Burette TITRONIC® 300 must be opened by authorised persons only; this means, for instance, that work on electrical equipment must only be performed by qualified specialists.

⚠ In the case of nonobservance of these provisions the Piston Burette TITRONIC® 300 may constitute a danger: electrical accidents of persons or fire hazard. Moreover, in the case of unauthorised intervention in the Piston Burette TITRONIC® 300 as well as in the case of negligently or deliberately caused damage, the warranty will become void.

Prior to switching the device on it has to be ensured that the operating voltage of the power supply and the Piston Burette TITRONIC® 300 matches the mains voltage. The operating voltage is indicated on the specification plate. Nonobservance of this provision may result in damage to the power supply Piston Burette TITRONIC® 300 or in personal injury or damage to property.

If it has to be assumed that safe operation is impossible, the Piston Burette TITRONIC® 300 has to be put out of operation and secured against inadvertent putting to operation. In this case please switch the Piston Burette TITRONIC® 300 off, pull plug of the mains cable out of the mains socket, and remove the Piston Burette TITRONIC® 300 from the place of work.

Examples for the assumption that a safe operation is no longer possible,

- the package is damaged,
- the Piston Burette TITRONIC® 300 shows visible damages,
- titrator Piston Burette TITRONIC® 300 does not function properly,
- liquid has penetrated into the casing.
- If the Piston Burette TITRONIC® 300 has been altered technologically or if unauthorized personnel tried or succeeded to open the instrument as attempt to repair it.

In case that the user operates such a device, all thereof resulting risks are on the user.

The Piston Burette TITRONIC® 300 must not be stored or operated in humid rooms.

For reasons of safety, the Piston Burette TITRONIC® 300 must only be used for the range of application described in the present operating instructions.

In the case of deviations from the intended proper use of the device, it is up to the user to evaluate the occurring risks.

⚠ The relevant regulations regarding the handling of the substances used have to be observed: The Decree on Hazardous Matters, the Chemicals Act, and the rules and information of the chemicals trade. It has to be ensured on the side of the user that the persons entrusted with the use of the Piston Burette TITRONIC® 300 are experts in the handling of substances used in the environment and in Piston Burette TITRONIC® 300 or that they are supervised by specialised persons, respectively.

⚠ During all work with titration solutions: Please wear protective glasses!

The Piston Burette TITRONIC® 300 is equipped with integrated circuits (EPROMs). X rays or other high energy radiation may penetrate through the device's casing and delete the program. For working with liquids, not being common titration solvents, especially the chemical resistance of the construction materials of the Piston Burette TITRONIC® 300 have to be considered (please also refer to chapter 1.1). For the use of liquids with high vapour pressure or (mixture of) substances not being mentioned in chapter 1.1 as allowed substances, the safe and proper operation of the Piston Burette TITRONIC® 300 has to be guaranteed by the user. When the piston moves upwards within the cylinder, a microfilm of dosing liquid or titration solution will always remain adhered to the inner wall of the cylinder, but this has no influence on the dosing accuracy. This small residue of liquid, however, may evaporate and thus penetrate into the zone underneath the piston, and if non-admitted liquids are being used, the materials of the Piston Burette TITRONIC® 300 may be dissolved or corroded (please refer also to chapter 8 "Maintenance and Care of the Piston Burette TITRONIC® 300").

2 Unpacking and First Operation

2.1 Unpacking and First Operation of the Piston Burette

The piston burette itself as well as all related accessory and peripheral parts have been carefully checked at the factory to ensure their correct function and size.

Please ensure that the small accessories are also removed in full from the packaging.

For the scope of delivery, please refer to the enclosed parts list.

The TITRONIC® 300 piston burette may be placed on any flat surface.

2.2 Connection and installing of Piston Burette and magnetic stirrer TM 50

The low voltage cable of the power supply TZ 1853 has to be plugged in to the 12 V socket „in“, (see Fig. 3 back panel, chapter. 2.3), on the back panel of the Piston Burette. Then plug the power supply into the plug socket.

Fig. 1a)

Place the power supply easily accessible in order to be able to remove the Piston Burette anytime easily from the power circuit.

The stirrer connects to the right side at the bottom and is locked in position by pushing it backwards (Fig. 1). In this way the power supply of the TM 50 stirrer is automatically established.

Fig. 1b)


Fig. 1b)

Fig. 2)

The stand rod TZ 1748 is screwed into the thread and the titration clamp Z 305 may now be mounted on the stand rod. (Fig.2). Instead of the magnetic stirrer TM 50, also the titration stand without stirring function (TZ3886) can be connected.

2.3 Connecting the piston burette – Combination with Accessories and Additional Devices

2.3.1 Back panel of the TITRONIC® 300 piston burette


Fig. 3

2.3.2 Connection ports of the TITRONIC® 300 piston burette

The burette TITRONIC® 300 offers the following interfaces:

Two RS232 interfaces, 4-channel (Mini-DIN):

- 1) For connection of a weighing balance and other devices from von SI Analytics (burettes, etc.)
- 2) RS1 for connection to the PC
- 3) One USB-A ("Host") interface for connecting USB devices such as a keyboard, printer, manual control unit, USB memory stick and USB hub
- 4) USB-B interface (mini USB) for connection to a PC
- 5) Connection of the external power supply TZ 1853
- 6) On/Off switch

2.3.3 Connecting a printer

Printers with a USB interface are to be connected to one of the two USB-A interfaces. These printers **have to** feature HP PCL emulation (3, 3 enhanced, 5, 5e, no GUI and GDI printers).

2.3.4 Connecting a USB device (manual controller, keyboard, memory device, hub)

The following USB devices can be connected to the USB-A interfaces:

- PC-keyboard
- TZ 3880 manual controller
- Printer
- USB storage devices, e.g. USB sticks
- USB hub
- USB barcode scanners

For connecting more than one USB device, a USB hub (e.g. TZ 3830) is required.

2.3.5 Connection of analytical balances

Analytical balances are to be connected to the RS232-2 port using an appropriate cable.

2.4 Setting the Language of the Country

The ex-factory default language setting is English. When the piston burette is switched on, the main menu will appear once the boot sequence is completed:


Fig. 4

Using <MODE>, followed by <System settings> you navigate to the system settings. The very first menu is to be used for setting the language of the country:


Fig. 5

The first menu is the language settings.


Fig. 6

Use <OK> to call the menu. Select the national language using the <↑↓> arrow keys, confirm it with <ENTER>/<OK>:


Fig. 7

The selected language will appear immediately. Pressing the <ESC> key twice will return the user to the main menu.

2.5 Dosing unit and Accessories



Fig. 8

- 1) TZ 2003 - drying tube
- 2) TZ 3282 - dosing hose without dosing tip and holding bracket
- 3) TZ 3802 - threaded cap with borehole GL 45, incl. adapter with 2 openings for drying tube and suction hose
- 4) TZ 1748 - stand rod
- 5) Z 305 - titration clamp
- 6) TZ 3620 - dosing hose with dosing tip and holding bracket: bracket = TZ 3875
- 7) TZ 3803 - 1 litre reagent bottle, brown
- 8) TZ 3656 - titration tip unit, blue
- 9) TZ 3801 - valve cover lid and TZ 3000 - 3/2-way valve
- 10) TZ 3130 - 20 ml dosing unit or TZ 3160 50 ml dosing unit
- 11) TZ 3283 - connection hose
- 12) TZ 3281 - suction hose

2.6 Installing the burette tip

The burette tip consists of the elements shaft with threaded clamping joint, hose and slip-on tip.

Burette tip - Sequence of assembly:

1. Cut of hose end evenly.
2. Slip parts of the threaded clamping joint on to the hose.
3. Guide hose through shaft.
4. Press the free hose end over the ring marking until it reaches the stop of the tip.
5. Push the tip with pressed in hose onto the shaft.
6. Hold tip firmly, and screw threaded clamping joint to the shaft.

2.7 Initial Filling or Rinsing of the Dosing unit

The dosing unit and the tubes are already mounted and ready-to-use. After the reagent bottle is connected, the initial filling of the dosing unit can be performed. While performing the initial filling or washing program, an adequately sized waste bin has to be placed below the titration tip.

On the main menu (fig. 9)


Fig. 9

Press <MODE> key and select Rinsing:


Fig. 10

Confirm the selection by pressing <OK>.

At this point you can select the number of rinsing cycles (Fig. 11). Initial filling requires a minimum of two rinsing cycles. You can stop the rinsing operation (Fig. 12 and 13) at any time by pressing <STOP> and then resume rinsing with <START>.


Fig. 11

The device fills first before the rinsing procedure starts:


Fig. 12

Fig. 13


Fig 14

Fig. 15

The rinsing program (Fig.12-15) can be canceled at all time with <STOP> and the continued with <START>. When the rinsing is finished, you can get back to the start menu by pushing 2 x<ESC>.

3 Working with the Piston burette TITRONIC® 300

3.1 Front Keyboard

Fig. 16

Apart from alphanumeric input (a-z, A-Z, 0-9) and a few other functions, almost all functions can be performed using the front keyboard.

<Mode>:	Methods selection, rinsing, system settings, balance data
<EDIT>:	Changing the current method, new method, copy and delete method, standard methods
<ESC>:	<ESC> will take you back to the previous menu level
<START/STOP>	Start and Stop of a current method
<FILL>:	Filling the burette
<↑>	Arrow-key up: Choosing of a menu or changing a number
<↓>	Arrow-key down: Choosing of a menu or changing a number
<→>	Arrow-key right: Changing position of the cursor, Enter

3.2 Display

The display consists of a graphical LCD display with a resolution of 320 x 240 pixels.


Fig. 17

3.3 Manual controller

The manual controller (Fig. 18) is needed for manual titration. It can also be used for starting dosage or other methods.

Fig. 18

Mode	Black key	Grey Key
Manual titration	Start of titration, single-step and continuous titration (please refer to chapter 3.6.1, manual titration)	Filling Stop of titration including evaluation
Dosage through Dosage method	Start dosage	Filling

3.4 External PC Keyboard (Optional)

Keys	Function
<ESC>	<ESC> will take the user to the previous level on the menu.
<F1>/<START>	Start of a selected method
<F2>/<STOP>	Stop of the current method
<F3>/<EDIT>	Change of the current method, new method, copy method, standard methods
<F4>/<FILL>	Fill the interchangeable unit
<F5>/ 	Display and modification of the balance data
<F6>/<MODE>	Selection of method, rinsing, system settings, balance data
<F7>/<SYS>	System settings (language selection, time/date ...)
<F8>/<CAL>	No function on the TITRONIC® 300
<F9>/+ / -	Change of sign
<F10>/<DOS>	Call dosing menu
Num/ Scroll and Lock/ Lock	Without function
Prt Sc and Sys Rq	Without function
<ESC>	Selection of the method-selection menu from the main menu. Elders: <ESC> will take you back to the previous level in the menu.
<↑> <↓> <←> <→>	Selection of individual menus and numeric values
0...9	Input of numeric values
<ENTER>	Confirmation of input parameters
<← Backspace >	Deletion of one input digit / an input character to the left of the blinking cursor
Letters, ASCII-symbols	Alphanumeric input possible. Uppercase and lowercase possible.
All other keys	Do not have any function

3.5 Menu Structure

There are 4 selection menus:

- Start or main menu
- Method parameters
- Method selection
- System settings

After power-up, the main menu is always the first menu to appear. The method displayed will always be the last method that was used (Fig. 19).


Fig. 19

Pressing <START> will result in the immediate execution of the method shown. <EDIT>/F3 will take you to the method parameters (Fig. 20).


Fig. 20

At this point you can

- modify the current method
- create a new method
- call and memorise standard methods
- copy or delete an existing method
- print a method (titration methods only)

Use the <↓> und <↑> keys to select the submenus, confirm your selection with <OK>/<ENTER>. <ESC> will take you back to the main menu.

<MODE>/F6 will take you to the method selection menu (Fig. 21).


Fig. 21

Existing methods can be selected by pressing the <↓> und <↑> keys and confirming the selection with <OK>/<ENTER>. Once the selection made, you will return to the main menu with the newly selected method. If no method is selected, <ESC> will also take you back to the main menu.

To navigate directly to the system settings (Fig. 22 and Fig. 23) you can use the <SYS>/F7 key; you can also navigate there through the method selection menu.


Fig. 22


Fig. 23

3.6 Main Menu

After power-up, the main menu is always the first menu to appear. The method displayed will always be the last method that was used (Fig. 24). In the present case it was a titration method.


Fig. 24

3.6.1 Manual Titration

Manual titration is always performed using the manual controller. Manual titration is impossible without the manual controller.


Fig. 25

Using <START> pressing the Black key on the manual controller will start the manual titration method.

Depending on the settings of the method, the system will prompt you for the sample description (Fig. 26) and the sample weight (Fig. 27). You can use an external PC keyboard (optional) to enter a 20-digit alphanumeric sample ID.


Fig. 26

Edit weight

003.37210g

 20 ml 07/20/14 19:39

Fig. 27

The balance data can be input using the front keyboard or the external keyboard. Please confirm your selection using <OK>/<ENTER>.

In the case of automatic takeover of the weighing-balance data, the sample weights will be read out of the balance data memory. If the memory does not contain any balance data, a message informing you of the absence of balance data will be displayed (Fig. 28).

Titration is running

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No balance data available. Wait for
automatic sample weight.

 20 ml 07/20/14 19:48

Fig. 28

Even at this moment, pressing the Print key on the balance will still cause the transfer of the balance data. After the input of the sample description and/or the sample weight/ sample volume, the following display will appear:

Titration is running

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0.000 ml

Speed 5 
 Stop
 20 ml 07/20/14 19:40

Fig. 29

You can control the metering rate with the black key of the "mouse". A single depression of the key will cause a step up to the first level. Depending on the size of the dosing unit, this corresponds to 0.0025 ml (20 ml) or 0.00625 ml (50 ml). Displayed are 0.003 or 0.006 ml from a single step.

If one keeps the black key depressed on the first level, titration will be continued at a low rate. If you press the black key fully down (2nd level) titration will proceed at a higher rate. The rate of the second level can be set in five stages using the <↓↑> arrow keys. These stages can also be changed during manual titration.

Titration is running
HCl titration

2.655 ml

Speed 3
Stop
20 ml

^ v |
STOP |
07/20/14 19:46

Fig. 30

Stage 5 corresponds to maximum titration speed. Speed is reduced by 50% each time.

Example: 20 dosing unit:

Stage 5	40.00 ml/min
Stage 4	20.00 ml/min
Stage 3	10.00 ml/min
Stage 2	5 ml/min
Stage 1	2.5 ml/min

As soon as manual titration is completed, press the <STOP/F2> key or approx. for 1 sec. the grey key of the "mouse". The titration result will be calculated and displayed and optionally printed on the connected printer or USB stick in PDF format:

End of titration
HCl titration

Consumption 1.803 ml
HCl 1.98 %

Back
Printing...!

ESC |
07/20/14 19:47

Fig. 31

<ESC> will take you back to the start menu way to start the next titration immediately. Filling of the interchangeable unit occurs automatically.

3.6.2 Dosage

To start a dosage method, please use the <START>/<F1> key or the black key of the manual controller ("mouse")

Main menu

0.000 ml
2.000 ml

dosing NaOH
Method parameter
Select method / system
20 ml

START |
EDIT |
MODE |
07/20/14 19:51

Fig. 32

Fig. 33


Fig. 34

The dosed volume will be briefly displayed before the display returns to the main menu.

The next dosage operation can be started immediately. Filling of the unit following dosage will not occur automatically, unless the maximum cylinder volume has been reached or the automatic filling option ("always") was activated. The unit can be filled at any time using <FILL>.

A dosing operation can also be performed without any dosing method with the <DOS>/<F10> key of the external keyboard:

Dosing volume

000.000 ml

Fig. 35

This is the point to input the volume which will be dosed following the confirmation with <OK>/<ENTER>.

Dose

1.201 ml

8.500 ml

20 ml STOP ESC 07/20/14 21:14

Fig. 36

Further dosages can be performed using <OK> or <ENTER.> Filling of the unit following dosage will not occur automatically here, unless the maximum cylinder volume has been reached. The unit can be filled at any time using <FILL>. <ESC> will take you back to the main menu.

4 Method Parameters

From the main menu shown in (Fig. 32), <EDIT> will take you to the method parameters:


Fig. 37

4.1 Method editing and new method

If you select <edit method> or <new method> you will be taken to the modification or new creation of a method. Selecting <new method> will always lead to the prompt for the input of a method name (Fig. 38). This prompt will not appear in the case of the modification of an already created method.


Fig. 38

The method name can contain up to 21 characters. Special characters are also possible. If no keyboard is connected, the method name being displayed has to be adopted (in the present case „Method 01“). Numbering of methods will occur automatically. Press <OK>/<ENTER> to confirm the input. The method name can be changed at any time. Please continue at this point with **Chapter 4.5**.

4.2 Default Methods

The <Default methods> item of the TITRONIC® 300 contains a series of ready-made standard methods which can be conveniently selected (Fig. 39).


Fig. 39

Once the selection made, you are directly prompted for the input of the method name (Fig. 40).


Fig. 40

The standard name may be adopted or modified. Subsequently, you will be taken to the <Change method parameters> item. Please continue at this point with **Chapter 4.5**.

4.3 Copy Methods

Methods can be copied or stored with a new name. If you select this function, the current method will be copied and you can include a new name (Fig. 41)


Fig. 41

A new name with the suffix [1] will be assigned automatically, so as to avoid the existence of 2 methods having the same name. Subsequently, you will be taken to <Change method parameters>. Please continue at this point with **Chapter 4.5**.

4.4 Delete Methods

If this function is selected, you will be prompted to know whether the current method is actually to be deleted. You have to reply **<Yes>** in explicit terms and also confirm this reply with **<OK>/<ENTER>**.


Fig. 42

4.5 Change Method Parameters

The input or modification of the method name was already described in **Chapters 4.1** and **4.3**.


Fig. 43

4.5.1 Method type

On the **<Method type>** you can select whether you wish to perform a manual titration or a dosage or whether you wish to prepare a solution (Fig. 44).


Fig. 44

The selection of the Method type will have an influence of the further parameterisation of the method. For instance, if you select the dosing mode, no selection of a formula will be available.

4.5.2 Result (only for manual Titration)

The <Result> menu offers the following possible settings:


Fig. 45

The **Result text** may contain up to 21 alphanumeric characters including special characters.


Fig. 46

Please confirm your input with <OK/<ENTER>.

4.5.2.1 Formulae for manual titration

On the **Formula Selection** submenu you can select the appropriate calculation formula:


Fig. 47

The following calculation formulae are available for manual titration:

Titration formula	Additional information
ml	Used to calculate the consumption in ml
$(ml-B) \cdot T \cdot M \cdot F1 / (W \cdot F2)$	Formula for calculating the concentration of a sample taking into account a blank value in terms of ml.
$(B-ml) \cdot T \cdot M \cdot F1 / (W \cdot F2)$	Formula for calculating the concentration of a sample taking into account a blank value in

	terms of ml. Reverse titration (examples. CSB, saponification number)
$(B \cdot F_3 - ml \cdot F_1) \cdot T \cdot M / (W \cdot F_2)$	Formula for calculating the concentration of a sample taking into account a blank value, including a multiplicative factor. Reverse titration.
$(W \cdot F_2) / (ml - B) \cdot M \cdot F_1$	Formula for calculating a titer (T) of a titration solution.

Legend of the abbreviations used:

ml:	Titration solution consumption in ml
B:	Blank value in ml, in most cases determined by titration
T:	Titer of the titration solution (e.g. 0.09986)
M:	Mol; Mol- or equivalent weight of the sample (e.g. NaCl 58.44)
F1	Factor 1, Conversion factor
F2	Factor 2, Conversion factor
W	Weight, sample weight in grams or sample volume in ml.

After selecting a formula, please confirm your selection with <OK>/<ENTER>:

Formula parameter
(ml-B)*T*M*F1/(W*F2)

B (Blank value)	0.0000 ml
T (Titre)	0.10000000
M (Mol)	36.46000
F1 (Factor 1)	1.0000
W (Amount)	auto
F2 (Factor 2)	1.0000

20 ml 07/20/14 22:25

Fig. 48

The values of the individual parameters of the selected calculation formula can now be input one by one.

Formula parameter
M (Mol)

00036.46000

20 ml 07/20/14 22:26

Fig. 49

4.5.2.2 Sample weight and volume (sample quantity)

Formula parameter
 $(ml-B)*T*M*F1/(V*F2)$

B (Blank value)	0.0000 ml
T (Titre)	0.10000000
M (Mol)	36.46000
F1 (Factor 1)	1.0000
W (Amount)	1.00000 g
F2 (Factor 2)	1.0000

20 ml 07/20/14 22:27

Fig. 50

Formula parameter
Amount

Weight manual

Weight automatic

Fixed weight

Manuel Volume

Fixed Volume

20 ml 07/20/14 22:28

Fig. 51

The Sample Quantity (W) item is used to select whether one is wishing to use a sample weight or a sample volume for titration or solution preparation.

You have the following options (Fig. 51):

- **Manual sample weight:** The sample weight is enquired by a prompt at the start of the method and manually input.
- **Automatic sample weight:** The sample weight is automatically transferred by a connected balance.
- **Fixed sample weight:** A fixed sample weight is input in g. This weight will then automatically be used for each start of the method without any sample weight being enquired by a prompt.
- **Manual sample volume:** The sample volume in ml is prompted at the start of the method and manually input.
- **Fixed sample volume:** A fixed sample volume is input in ml. this volume will then automatically be used for each test of the method without any sample volume enquired by a prompt.

4.5.2.3 Formula unit

The formula unit can be selected in the **Unit** submenu.

Unit 1
HCl titration

None

ml

%

ppm

g

mg

20 ml 07/20/14 22:40

Fig. 52

Once the selection made (e.g. %), the unit will also be displayed as piece of information on the display.


Fig. 53

4.5.2.4 Decimal digits

To conclude, it is possible to determine the number of decimal digits from 2-6. The standard setting is 2.

4.5.2.5 Global Memories

If a titration result should be used again later, e.g. the factor or titer of a solution or a blank value, it can be saved automatically. The creation of a global memory is possible when using an external keyboard only. The creation of a global memory can be managed with the system settings or by typing Shift + F5 on the external keyboard. This leads to the <Global memory>:


Fig. 54

A global memory can be added by pressing F3:


Fig. 55


Fig. 56

MO1 can either be accepted or substituted by a designation of a blank value or a titer.


Fig. 57

This value can be used with the formula.


Fig. 58

The blank value, being titrated previously is always taken into account automatically.


Fig. 59

4.5.3 Dosing parameters

Edit method parameter
HCl titration

Method name
Method type man
Result

Dosing parameter

Sample ID man
Documentation GLP

^ V | OK | ESC |

20 ml 07/21/14 12:40

Fig. 60

The dosing parameters (dosing speed, filling speed and max. dosing/titration volume) are determined for each method. This applies to both types of methods such as manual titration and dosage.

Edit dosing parameter
HCl titration

Dosing speed 100.00 %

Dosing speed 40.0 ml/min
Filling speed 20 s
Max. titration volume 50.000 ml

^ V | OK | ESC |

20 ml 07/20/14 22:50

Fig. 61

The dosing speed can be set as a function of the interchangeable unit in terms of ml/min from 0.1 to 100 ml/min.

Dosing unit	Max. dosing speed [ml/min]
20 ml	40
50 ml	100

The filling speed can be set in terms of seconds from 20 to 999. The standard setting of this value is 30 seconds. For diluted aqueous solutions the filling speed can be six to 20 seconds. For non-aqueous solutions the filling speed should be set to the 30 seconds. In the case of highly viscous solutions such as concentrated sulphuric acid the filling speed should be further reduced down to 40 - 60 seconds.

Depending on the method type, the (maximum) the living volume or titration volume can be set to 999.999 or even 9999.999.

4.5.4 Sample identification

In the manual titration it is possible to input a sample identification. The possible input includes manual, automatic or no sample identification at all.


Fig. 62

For the manual sample ID, a prompt for the sample ID will always be displayed at the start of the method (Cp. also chapter 3.6, Main menu). For an 'automatic' sample ID there will be selected a master ID (in the current case this is water, cp. Fig. 68), which will then automatically be numbered starting on 01.


Fig. 63

After a new power-up, numbering will resume with 01.

4.5.5 Documentation


Fig. 64

Two different format settings are available for documentation on a printer or USB memory stick: short and GLP:


Fig. 65

Method type	Short documentation	Standard documentation	GLP-Documentation
Manual titration	Method name, date, time, sample description, sample weight/sample volume, results and calculation formula	N/A	Same as 'Short documentation' + plus method contents
Dosing	Only method printout, Method name, date, time, dosing parameters	N/A	N/A

The documentation can also be stored in PDF and CSV-files on a connected USB stick.

5 System Settings

Main menu

0.000 ml
0.000 ml

Methode 01 START
 Method parameter EDIT
 Select method / system MODE
 20 ml 07/21/14 12:47

Fig. 66

From the main menu (Fig. 66), press <MODE> followed by <System settings> (or with the external keyboard connected, press F7) you get to the system settings:

System settings

Language settings
 Reagents / dosing unit
 Global memory
 RS232 Settings
 Printer PDF
 Date/time
 Reset ▼
^ v OK ESC
 20 ml 07/21/14 12:47

Fig. 67

Setting the national language was already described in **Chapter 2.5**.

5.1 Reagents/Dosing Unit

In this menu you can set the size of the dosing unit (20 or 50 ml), perform an exchange of the dosing unit and enter reagents data, being stored during manual titration in the GLP documentation.

System settings
 Reagent
Unit size 20 ml
 Reagent
 Concentration 1.00000
 Conc. determined at --
 Expire date 07/11/14
 Opened/compounded 07/11/14 ▼
^ v OK ESC
 20 ml 07/21/14 12:52

Fig. 68

5.1.1 Dosing unit exchange

Confirm with <OK/ENTER> the size of the dosing unit.


Fig. 69

Chose <Dosing unit exchange>.

Warning: The exchange of the dosing unit starts without any warning. Therefor please place the titration tip in a reagents bottle or beaker.

The piston is moving up to 85%.

System settings
Dosing unit exchange

Dosing unit moves up

Fig. 70

Then the system prompts to unlock the dosing unit.

System settings
Dosing unit exchange

Please unlock dosing unit

Fig. 70

Then please unlock the dosing unit as shown in fig. 71:


Fig. 71:


Fig. 75

After this procedure the dosing unit is moving back down again. Please lock the dosing unit.



Abb. 76

The following reagents data can be chosen:

- Unit size 20 or 50 ml (adjustable)
- Reagent name: (default: blank)
- Concentration: (default: 1.000000)
- Concentration determined on: (default: actual date)
- To be used until: (default: actual date)
- Opened/Produced on: (default: actual date)
- Test according to ISO 8655: (default: actual date)
- Charge description: (default: blank)
- Last modification: (default: actual date)


Fig. 77

5.2 Global Memories

Please refer to chapter 4.5.2.5.

5.3 RS232 Settings

The <RS232 settings> item can be used to determine the device address of the TITRONIC® 300 and set the parameters of the two RS232 interfaces separately:


Fig. 78

The device address can be set from 0 – 15. Address 1 is the default setting:


Fig. 79


Fig. 80

The baud rate is preset to 4800. It may be set to 1200 – 19200:


Fig. 81

The parity setting can be selected amongst <No>, <Even> and <Odd>. <No> is the standard setting.


Fig. 82

You may select between 7 and 8 data bits. 8 bits is the standard setting.


Fig. 83

The stop bits can be set to 1, 1,5 and 2. The default setting is 1.


Abb. 84

Selecting <Reset RS Parameters> will reset the RS232 parameters to the factory settings.

5.4 Connection of Printers

The results and methods can be printed on the following media:

- HP PCL compatible printer (A4), colour (chromatic) and monochrome (e.g. laser printer)
- Seiko DPU S445 (Thermo paper 112 mm width)
- On the USB stick in PDF format

To connect the printers to the burette please use the USB socket. When printing, please check whether the correct printer is connected. It is not possible to print "HP" printer layouts on a Seiko thermal printer or vice versa. The printer settings should always be checked and adjusted after changing the printer.


Fig. 85

Only one printer at a time can be connected, because automatic printer recognition is not supported. HP PCL is the default setting. If you select "Print PDF", please make sure that a USB stick is connected to the device. Print PDF is the default setting.

5.5 Date and Time

The factory time setting is Central European Time. This setting may be changed, where necessary:


Fig. 86

5.6 RESET

RESET will reset all settings to the factory setting.

Please note: All methods will also be deleted. So please print the methods or export/copy them to a **connected USB storage medium (refer to chap. 5.9)**.

The RESET has to be confirmed separately once again:


Fig. 87

5.7 Device Information

<Device Information> contains information about

- the current software version
- the serial number of the device
- printer driver and update version
- device address set
- number of measurements (Starts of a method)
- a number of strokes/filling cycles


Fig. 88

Please hold this device information ready for service purposes.

5.8 System Sounds

The system tone can be switched on and off.


Fig. 89

5.9 Data exchange

All methods with all parameter settings and global memories can be stored and restored on a connected USB-memory. It is also possible to transfer the settings from one titrator to another one. The backup will be started with **Settings backup**:


Abb. 90

Back settings is displayed during the backup in blue:


Abb. 91

After a Reset or a maintenance case it is possible to restore the backup with **Restores settings**:

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5.10 Software Update



Fig. 94

An update of the device software requires a USB stick containing a new version. For this operation, the two files that are needed have to be located in the root directory of the USB stick:

Plug the USB stick into a free USB-A port, wait for some seconds, and then select the Software Update function. The valid software updates will be shown on the display. In the present case this is Version „14_29“ from week 29, year 2014.


Fig. 95

After starting the update using <OK/ENTER>, next thing to appear is the following graphic:


Fig. 96

which will change after a few seconds to the following display:


Fig. 97


Fig. 98

Upon completion of the update (approx. 2-3 minutes), the device will shut down the software completely and proceed to a new start.

Please note: During the update, the device must not be switched off.

Important: In the course of an update, the methods will not be deleted! You can continue to use them. If no valid update file is stored on the USB stick, the following message will appear:

Software Update

Software version: 1429

No update found

Fig. 99

6 Connection of Analytical Balances

As it often happens that the sample is weighed in on an analytical balance, it makes sense to connect this balance to the TITRONIC® 300. To connect the balance to the TITRONIC® 300, the balance must have a RS-232-C-interface and the connection cable must be configured accordingly. For the following types of balances there are already assembled connection cables:

Balance	TZ-Number
Sartorius (all types with 25 pin RS)	TZ 3092
Mettler, AB-S, AG, PG	TZ 3099
New Sartorius balances with USB-Interface (via RS-converter)	TZ 3099
Precisa XT-Series	TZ 3183
Kern with 9-pole RS232	TZ 3180

For all other types of balances it is possible to obtain an already assembled connection cable (on demand). For this we need detailed information about the RS-232-C-interface of the balance used.

The connection cable is to be connected to the RS-232-C-interface 2 of the TITRONIC® 300. This side of the connection cables always consists of a 4-pole mini-plug. The other side of the cable can, depending on the type of balance, be a 25-pole plug (Sartorius) or a 9-pole plug (Mettler AB-S).

In order to allow the balance data to be sent to the TITRONIC® 300, the data transmission parameters of the titrator and the balance must correspond to each other. Additionally, it is necessary to carry out some more standard settings on the side of the balances:

- The balance is to send the balance data via RS-232-C only by means of a print command.
- The balance is to send the balance data only after the display standstill.
- The balance should never be set to 'automatic sending' and/or 'send continuously'.
- 'Handshake' on the balance must be set to 'off', or even 'Software Handshake' or 'Pause'.
- No special characters such as **S** or **St** are allowed to be used as prefix in the balance data of the balance data string. In such a case it might be possible that the TITRONIC® 300 cannot process the balance data correctly.

After you have connected the balance with the appropriate cable to the TITRONIC® 300 and have adjusted all settings in the balance software, and possibly in the TITRONIC® 300, you can now test the data transfer of the balance very easily. Start the one method. Confirm the sample designation. Then, the display asks you:

- a) To press the print-button at the balance → Parameters to 'weighted sample automatically'
- b) To enter the weighted sample → then the parameters are still set to 'weighted sample manually'

Put an object onto the balance and press the print button. After the standstill of the balance display there will be beep at the TITRONIC® 300 and the transmitted balance data appear:

- After approx. 5 sec. in the display and the display changes automatically into the measuring display.
- The weighted sample must again be confirmed with <Enter> or <F1>.

6.1 Balance data editor

Pressing the die **<F5/balance symbol >** function key will invoke the so-called balance data editor. A list with the existing balance data will appear:


Fig.100

The balance data can be edited one by one. Following a change, a star will appear opposite the weighed-in quantity:


Fig. 101

Weights may be deleted or added individually. It is also possible to delete all weights at one stroke.


Fig. 102

If no balance data is available, the “No balance data found” message will appear.:


Fig. 103

7 Data Communication via RS-232- and USB-B interface

7.1 General Information

The burette TITRONIC® 300 has two serial RS-232-C interfaces to communicate data with other devices. By means of these two interfaces it is possible to operate several devices on one computer (PC) interface.

In addition to that, the TITRONIC® 300 also has an alternatively USB-B interface, which can only be used to connect a PC.

RS-232-C-1 establishes the connection to a connected computer or to the previous device of the “Daisy Chain”. At the RS-232-C-2 it is possible to connect additional devices (Daisy Chain Concept).

PIN assignment of the RS-232-C interfaces:	PIN-No.	Meaning / Description
	1	T x D Data output
	2	R x D Data input
	3	Digital mass

7.2 Chaining multiple devices — „Daisy Chain Concept“

In order to activate several devices in a chain individually, each device must have an own device address. For this it is at first necessary to establish a connection from the computer to the RS-232-C interface 1 of the first device in the chain by means of a RS-232-C data cable, e.g. Type No. TZ 3097. With the additional RS-232-C data cable, Type No. TZ 3094, the RS-232-C- interface 2 of the first device is connected with the RS-232-C-interface 1 of the second device. At interface 2 of the second device it is possible to connect an additional device.

The TITRONIC® 300 can also be connected via USB cable TZ 3480 (type A (M) – type B (M), 1.8m). It is also possible to connect the TITRONIC® 300 via a mini USB to a USB interface of a PC. To accomplish this connection, a driver has to be installed on the PC. Then the USB-B interface takes over the function of the RS232-1 interface. Please contact SI Analytics for acquiring the software.

The address always consists of two characters: e.g. address 1 of the two ASCII- characters <0> and <1>. The addresses can be set from **00** to **15**, i.e. 16 possibilities. It must be ensured that the devices in a chain have different addresses. If a device is addressed with its address, this device will process this command without sending it to another device. The reply to the computer has also an own address. The addresses are allocated as described in  **Chapter 5.2**.

The burette TITRONIC® 300 receives commands from a PC at the interface 1 (USB- B) if the computer knows the address. It also sends the answer via this interface. If the address of the incoming command does not match the device address, the complete command will be forwarded to interface 2. Interface 2 is connected to interface 1 of another device. This device checks the address as well and reacts to the command as the first TITRONIC® 300 did before.

All information (data strings) which arrive at interface 2 of the burette TITRONIC® 300 will immediately be send to the computer via interface 1 (or USB-B interface). Thus, the computer receives the data of all devices. In practice it is possible to connect up to 16 devices to one computer- (PC-) interface.

7.3 Instruction Set for RS-Communication

The commands consist of three parts: Address two-digit aa, e.g.: 01
 Command e.g.: DA
 Variable, if necessary e.g.: 14
 and end of command <CR> <LF>

Every command must be completed with the ASCII - sign <CR> and <LF> (Carriage Return and Line Feed). Only if the respective action has ended the answers will be returned to the computer.

Example: The command to dose 12.5 ml shall be sent to the burette TITRONIC® 300 with the address 2.

The command consists of the characters: 02DA12.5<CR LF>

In detail: 02 = Device address
 DA = Dosage command with filling and zero points of the display
 12.5 = Volume in ml to be dosed
 <CR LF> = Control character as command end

Command	Description	Reply
aaAA	automatic allocation of device address	aaY
aaMC1...XX	choosing a method	aaY
aaBF	„filling burette“. Aufsatz wird gefüllt.	aaY
aaBV	output of dosed volume in ml	aa0.200
aaDA	dose volume without filling, with adding the volume	aaY
aaDB	dose volume without filling, reset of the volume	aaY
aaDO	dose volume with filling, without adding the volume	aaY
aaGDM	dosing speed in ml/min	aaY
aaGF	filling time in seconds (min is 20, default 30)	aaY
aaES	„ESC“ function one step backwards	aaY
aaEX	„EXIT“ function.back to main menu	aaY
aaGDM	dosing speed in ml/min (0.01 – 100 ml/min)	aaY
aaGF	filling time in sec (adjustable 20 – 999 seconds)	aaY
aaGS	output serial no. Of device	aaGS08154711
aaLR	output report (short report)	aaY
aaLI	output method content	
aaLO	output documentation (as configured)	
aaRH	request of identification	aaIdent: TITRONIC 300
aaRC	send last command	aa“last command“
aaRS	report status possible answers are: „STATUS:READY“ for ready „STATUS:dosing“ dosing „STATUS:filling“ filling „ERROR:busy“ if no interchangeable unit has been attached	aaStatus:“text“
aaSM	start selected method	aaY
aaSEEPROM	EEPROM reset to factory defaults	aaY
aaSR	stop the actual function	aaY
aaVE	Version number of the software	aaVersion:

8 Maintenance and Care of the TITRONIC® 300 Piston Burette

The preservation of the proper functioning of the piston burette requires testing and maintenance work to be performed on a regular basis.

Regular inspections are essential prerequisites for the correctness of the volume and the proper functioning of the piston burette.

The accuracy of the volume is determined by all chemicals-carrying components (piston, cylinder, valve, titration tip and hoses). These parts are subject to wear and tear, i.e. they are or wearing parts, respectively. The piston and cylinder are subject to particular strain, hence they require special attention.

Heavy strain:

Use of e.g. concentrated solutions, reagents and chemicals (> 0,5 mol/L); chemicals attacking glass, such as fluorides, phosphates, alkali solutions; solutions with a tendency to crystallising out; Fe (III) chloride solutions; oxidising and corroding solutions such as iodine, potassium permanganate, Cer (III), Karl-Fischer titration agent, HCl; solutions with a viscosity of > 5 mm²/s; frequent, or even daily use.

Normal strain:

Use of solutions, reagents and chemicals (up to 0.5 mol/l) which do not attack glass, crystallise out or corrode.

Interrupted use :

If the dosing system is not in use for more than two weeks, we recommend emptying and cleaning the glass cylinder and all hoses [6]. This applies in particular under the operating conditions referred to in the "Heavy strain" section. If this recommendation is not adhered to, the piston of the valve may become leaking, this may result in damage to the piston burette.

If the liquid is left within the system, you will also have to reckon with corrosion and an alteration of the solutions used over time, which includes e.g. crystallisation. Considering that as of the state of the art there are no plastic hoses available for the use in titration equipment which would be perfectly free of diffusion phenomena, particular attention is to be paid to the range of the hose lines.

We recommend the following inspection and maintenance work	Heavy strain	Normal strain
Simple cleaning: <input type="checkbox"/> Wiping off splashed chemicals from the outer surface. [1]	Whenever required in operation	Whenever required in operation
Sight check: <input type="checkbox"/> Check for leakage in the area of the dosing system. [2] <input type="checkbox"/> Is the piston tight? [3] <input type="checkbox"/> Is the valve tight? [4] <input type="checkbox"/> Titration to clear? [5]	Weekly, when putting back into operation	Monthly, when putting back into operation
Basic cleaning of the dosing system: <input type="checkbox"/> All parts of the dosing system to be cleaned separately. [6]	Every three months	Whenever necessary
Technical inspection: <input type="checkbox"/> Check for air bubbles in the dosing system. [7] <input type="checkbox"/> Visual inspection <input type="checkbox"/> Check of the electrical connections. [8]	Semi-annually when putting back into operation	Semi-annually when putting back into operation
Verification of the volume according to ISO 8655: <input type="checkbox"/> Perform basic cleaning <input type="checkbox"/> Inspection according to ISO 8655 Part 6 or Part 7. [9]	Semi-annually	Annually

Please note: Depending on the respective application, there may be different specifications for the entirety of the inspection and maintenance work to be performed. The individual intervals may be extended if no complaints occur, but they will have to be shortened again as soon as any problem has arisen.

The inspection of the metrological reliability including maintenance work is offered as a service (including a manufacturer's certificate, if so ordered). In this case the titration device is to be sent in. Please contact the service (see backside of this manual).

Detailed description of the inspection and maintenance work:

- [1] Wipe off using a soft cloth (and some water with a normal household detergent).
- [2] Leaking connections can be identified by moisture or crystals at the threaded connections of the hoses, at the sealing lips of the piston inside the dosing cylinder or at the valve.
- [3] If any liquid becomes visible below the first sealing lip, it has to be checked at short timely intervals whether any liquid will build up under the second sealing lip, too. In this case both the piston and the glass cylinder have to be replaced immediately. It is easily possible that in operation small liquid droplets build up under the first sealing lip, but they may also disappear again. This phenomenon alone is no reason for replacement.
- [4] The valve has to be removed from its housing for inspection. In this process, the hoses remain connected to the valve. Please check for moisture underneath the valve. When reinserting the valve, please make sure that the small cam at the rotating axis is fitted into the corresponding groove again.
- [5] The titration tip must be free of sedimentation or crystals which might obstruct the dosing process or falsify the results.
- [6] Remove the cylinder, take the valve out of the valve housing, unscrew the hoses and then rinse all parts carefully with distilled water. For the assembly of the cylinder, hoses and other parts of the interchangeable unit, please refer to the operating instructions.
- [7] Dose one burette volume, then refill. Air bubbles will gather at the tip of the cylinder and in the titration hose where they can be detected easily. If bubbles become visible, please re-tighten all connections finger tight, and then repeat dosing. If air bubbles still remain within the system, [6] please check the valve and replace the hose connections. The air bubbles may also occur at the interface between the sealing lip of the piston and the cylinder. If a reduction of the filling speed will not do, the dosing unit has to be replaced.
- [8] Check the electrical plug contacts for corrosion and mechanical damage. Defective parts have to be repaired or replaced by new parts.
- [9] Please refer to the application "Burette inspection according to ISO 8655 Part 6".

9 Storage and transportation

If the TITRONIC® 300 Piston Burette or the interchangeable units have to be stored over some time, or to be dislocated, the use of the original packing will be the best protection of the devices. However, in many cases this packing will not be available anymore, so that one will have to compose an equivalent packaging system. Sealing the lower section in a foil is hereby recommended.

The devices should be stored in a room with a temperature between +10 and +40°C, and the (relative) humidity of the air should not exceed 70 %.

If the interchangeable have to be stored over some time, or to be dislocated, the fluids inside the system, especially aggressive solution have to be removed (please refer also to chapter 8. „Maintenance and Care of the TITRONIC® 300 Piston Burette“).

10 Recycling and Disposal

Please observe the applicable local or national regulations concerning the disposal of "waste electrical and electronic equipment".

The present piston burette and his packaging are manufactured as far as possible from materials which can be disposed of environmental-friendly and recycled in a technically appropriate manner. If you have any question regarding disposal, please contact the service (see backside of this manual).

i The main printed board carries a lithium battery. Batteries should not to be disposed of with the normal domestic waste. They will be taken back and recycled or disposed of properly by the manufacturer at no cost.

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Xylem Analytics Germany Sales GmbH & Co.KG