Lovibond® Colour Measurement

Tintometer® Group





Lovibond® PFXi-195 Series

Dark Oils & Fats • Liquid Chemicals & Industrial Oils
Petroleum Oils & Fuels • Beers, Malts & Caramels • Pharmaceutical Solutions
Industrial Oil & Surfactants • Sugar Solutions, Syrups & Honey
Optically Clear Samples • Water & Waste Waters

Accurate, efficient spectrophotometric colorimeters

www.lovibondcolour.com

Lovibond® PFXi-195 Series

- RCMSi Remote Calibration & Maintenance Service via internet
- USB & LAN Interface
- Large clear LCD screen for ease of use and ability to display graphs and data in a wide variety of languages
- Extensive yet flexible choice of standard colour scales
- Remote upgrade facility for adding scales once in service
- Allows calculation and description of off-hue status
- Gives closest match to stored references
- Generates a customised colour scale from reference samples
- Robust steel construction with excellent chemical resistance
- Easy maintenance with removable sample chamber
- Supplied with colour control software for data analysis
- Output conforming to GLP including date, time, sample & user ID
- Accommodates a range of sample cells and tubes

Objective Colour Data for Process Control

The Lovibond® PFXi-195 Series is a highly efficient colorimeter for optically clear samples that meets the demand for consistent and reliable colour data. It removes all subjectivity involved in colour grading, supplying unbiased readings that are unaffected by operator or environment.

Comprehensive Selection of Standard Colour Values

The PFXi-195 Series of colorimeters automatically measure colour and display the results directly, either according to the traditional grading scales that have been widely adopted as industry standards for colour assessment and control, or in terms of internationally recognised CIE values and spectral data.



Colour Scale	References	Scope	Range
ADMI (full spectrum and tristimulus filter)	American Standard Methods 2120 E	Coloured waters and tinted liquids	1 – 500
Acid Wash Color	ASTM D848	Quality testing of industrial aromatic hydrocarbons	1 – 14
ASBC		American standard for colour grading of beers; derived from EBC Colour	1.2 – 10.6 (extended range by dilution and reduction in path lengtl
ASTM Color	ASTM D1500, D6045, ISO 2049	A wide range of petroleum products including lubricating oils, heating oils and diesel fuel oils	0.5 – 8 units
Chinese Pharmacopoeia Colour Series ¹⁾	CP Appendix IX A	Pharmaceutical solutions	YG1 – 10; Y1 – 10; OY – 10; OR1 – 10; BR1 – 10
EBC Colour	Analytica	Beers, malts and caramels and similarly coloured liquids. Based either on absorption at 430 nm or CIE x y chromaticity co-ordinates	2 – 27 units (extended range by dilution and reduction in path length
European Pharmacopoeia Colour Series	Ph. Eur. Method 2.2.2	Pharmaceutical solutions	R1 – 7; Y1 – 7; B1 – 9; BY1 – 7; GY1 – 7
FAC Color	AOCS Cc 13a-43	Approved by the Fats Analysis Committee of the American Oil Chemists Society for grading dark coloured oils, fats and tallows.	1 – 45 (odd numbers)
Gardner Color	ASTM D1544, D6166, AOCS Td 1a,	Chemicals and oils ranging from pale yellow to red, such as resins, varnishes, drying oils, lecithins and fatty acids	1 – 18 units
Hess-Ives Color Units	DGK F050.2	Chemicals and surfactant liquids	
Honey Colour (Pfund Equivalents)		Commercial honeys, ranging from pale yellow through amber to deep red	0 – 115 mm Pfund
ICUMSA Colour	ICUMSA GS1-7, ICUMSA GS2/3-9	Sugar solutions & syrups	
Iodine Colour	DIN 6162	Solvents, plasticisers, resins, oils and fatty acids ranging from yellow to brown	1 – 500 units
Klett Color (blue filter KS-42)	AOCS Dd 5-92	Detergents and surfactants	0 – 1000 units
Pt-Co/Hazen/APHA Scale	ASTM D1209, D5386	Water and other clear liquids such as plasticisers, solvents and petroleum spirits	0 – 500 mg Pt/l
Rosin, US Naval Stores ¹⁾	ASTM D509	Rosins varying in colour from yellow to reddish orange	XC – D + FF
Saybolt Color	ASTM D156, D6045, JIS K 2580	Light coloured petroleum products including aviation fuels, kerosine, white mineral oils, hydrocarbon solvents and petroleum waxes	–16 (darkest) to +30 (lightest)
Series 52 (Brown)		Beers, whiskies and sugar solutions	1– 38 units
Yellowness Index ¹⁾	ASTM D1925, E 313	Determination of the degree of yellowness under daylight illumina Calculated from X Y Z tristimulus values	ition.
US Pharmacopoeia Color	USP (631) Color and Achromicity	Pharmaceutical solutions	A – T
CIE Values	ASTM E308	$X\ Y\ Z$ tristimulus values; $x\ y\ Y$ chromaticity co-ordinates; CIE $L^*a^*b^*$ colourspace; ΔE colour difference; L^*C^*h colour space ²⁾ , Hunter $L\ a\ b$ colour sp	Defined by spectrum locus bace ²⁾
Spectral data (420 – 710 nr	m)	Transmittance (full spectrum and at specified wavelengths) Optical density (full spectrum and at specified wavelengths)	0 – 100% 0 – 2.5

included as standard on PFXi-195/9 only

Version	Application	Standard Colour Scales
PFXi-195/1	Liquid chemicals & industrial oils	Pt-Co/Hazen/APHA, Gardner, Iodine, CIE values, spectral data
PFX <i>i</i> -195/2	Petroleum oils & fuels	Saybolt, ASTM Color, Pt-Co/Hazen/APHA, CIE values, spectral data
PFX <i>i</i> -195/3	Dark oils & fats	FAC, Gardner, CIE values, spectral data
PFX <i>i</i> -195/4	Beers, malts & caramels	EBC (CIE & 430 nm), ASBC (CIE & 430 nm), Series 52 (Brown), CIE values, spectral data
PFXi-195/5	Pharmaceutical solutions	European Pharmacopoeia (includes automatic selection of colour series), US Pharmacopoeia, Pt-Co/Hazen/APHA, CIE values, spectral data
PFXi-195/6	Industrial oils & surfactants	Klett Color (blue filter KS-42), Pt-Co/Hazen/APHA, Hess-Ives Color Units, CIE values, spectral data
PFXi-195/7	Sugar solutions, syrups & honey	ICUMSA Colour (420, 560, 710 nm), Honey Colour (Pfund Equivalents), Series 52 (Brown), CIE values, spectral data
PFXi-195/8	Waters & waste waters	ADMI (spectral & tristimulus filter methods), Pt-Co/Hazen/APHA, CIE values, spectral data
PFXi-195/9	Transparent samples	CIE values (standard selection), L*C*h colour space, Hunter L a b colour space, spectral data

Colour scale upgrades give the flexibility to meet individual requirements, enabling additional colour scales to be added to standard instrument versions either at the time of order or remotely once the instrument is in service. For product types that are incompatible with standard colour scales, the PFX*i*-195 allows users to build up a customised scale from a series of reference samples and then to obtain a closest match to the stored references.

RCMSi Technology (Remote Calibration & Maintenance Service via internet)

The PFXi range includes new RCMSi technology (Remote Calibration & Maintenance Service via internet). This unique feature allows a calibrated measurement to be taken using an ISO 17025 certified liquid standard. The spectral response is transmitted to the Lovibond® secure calibration server. On completion of this procedure, a traceable calibration certificate is made electronically available to the user. The RCMSi system not only gives the user added confidence that the instrument is working within tight ISO standards but also reduces the need for expensive on-site servicing and preventative maintenance of the unit.

Colour Testing Made Simple

The PFXi-195 Series provides easy to use, automatic instruments that require no special skills to operate. The built-in menu guides users through the selection of operating parameters such as colour scale. Thereafter, readings are made with a single key press, taking less than 20 seconds to complete. Built-in help files are available.

Easily Customised to User Specifications

Adaptable software and design allow users to configure the PFX*i*-195 Series to their requirements. Operators can set the language for display, program the PFX*i*-195 to show only those scales of interest or restrict access to the menu system. As well as standard colorimetry cells, the PFX*i*-195 can be used with a range of tubes and standard, flow-through and disposable spectrophotometer cells.



Calculation and Description of Off-Hue Status

The PFX*i*-195 off-hue status is a useful facility that reveals whether the sample colour is characteristic of the selected scale. Where appropriate, it includes a description of hue difference (eg. redder, greener), relative saturation (stronger or weaker) and an off-hue factor (a relative measure of the distance away in colour space of the sample colour from the 'true' colour scale).

Suited to Laboratory or Production Environments

Comprehensive facilities for colour management make the PFX*i*-195 Series an ideal choice for the laboratory. However, with excellent calibration stability, password protection for tamper proof control and simple operation, the PFX*i*-195 also supports the migration of quality control to the manufacturing area, making it a cost effective option for dedicated production testing. For easy maintenance, the Lovibond® PFX*i*-195 Series includes a robust steel sample chamber, which is simply removed and cleaned if a spillage occurs, and the precision filament lamp is easily accessed and changed from outside the instrument.

Optimised Use of Colour Data

Data sets can be saved in the instrument, printed out or automatically downloaded to a PC where they can be processed and stored for future analysis, traceability and trend monitoring. The colour control software supplied with the PFX*i*-195 enables the generation of spectral and CIE diagrams as well as analysis of spectral data. It also permits direct control of the instrument from the computer.

Confidence in Colour Measurement

For regular conformance checking, certified colour reference standards are also available for routine calibration and verification of test data. Conformance filters are supplied with a Certificate of Conformity stating their colour values and confirming that they have been manufactured and inspected under the control of The Tintometer Ltd's ISO 9001 quality system. Certified colour reference standards are supplied with full traceability to internationally recognised standards, either ISO/IEC 17025 (ASTM Color, Saybolt and Gardner Color) or the ISO 9001 quality system (Pt-Co Units).

Lovibond® PFXi-195

PFXi-195 Technical Specification

Measuring principle	9 interference filters
Spectral response	420 – 710 nm
Bandwidth	10 nm
Repeatability	on measurements of de-ionised water
Chromaticity (x y)	+/- 0.0004
ΔΕ	0.4
Measurement time	Less than 30 seconds
Baseline calibration	Single key press; fully automated
Light source	5 Volt, 10 Watt Tungsten Halogen
Illuminants	CIE Illuminant A, B, C, D65
Observer	2°, 10°
Path length	0.1 – 50 mm
Interfaces	USB, LAN, RS232
Data storage	100,000+ Measurements
Input voltage	Universal (90 – 240 Vac), via external power supply
Compliance	CE, RoHS, WEEE
Display	240 x 128, back-lit graphic display, (white on blue)
Keypad	23 key tactile membrane keypad;
	washable polyester with audible feedback
Languages	English, French, German, Italian, Portuguese,
	Russian, Spanish, Chinese and Japanese
Instrument housing	Fabricated steel with tough, textured paint finish
Dimensions	Width 435, depth 195, height 170 (mm)
Weight	6.8 kg

Accessories

Version	Optical Glass Cells
PFXi-195/1	1 x 10 mm, 1 x 50 mm
PFXi-195/2	1 x 33 mm, 1 x 50 mm
PFXi-195/3	1 x 10 mm
PFXi-195/4	1 x 10 mm, 1 x 25 mm
PFXi-195/5	1 x 40 mm, 1 x 50 mm
PFXi-195/6	1 x 10 mm, 1 x 40 mm, 1 x 50 mm
PFXi-195/7	1 x 10 mm, 1 x 50 mm
PFXi-195/8	1 x 10 mm, 1 x 50 mm
PFXi-195/9	1 x 10 mm

^{*} Whenever cells are subjected to thermal shock, only borosilicate glass cells should be used.

Ordering Information Order Code PFXi-195 Instruments

PFXi-195/1 (Pt-Co, Gardner, Iodine)	1371951
PFXi-195/2 (Saybolt, ASTM Color, Pt-Co)	1371952
PFXi-195/3 (FAC, Gardner)	1371953
PFX <i>i</i> -195/4 (EBC (CIE & 430 nm), ASBC (CIE & 430 nm), Series 52)	1371954
PFX <i>i</i> -195/5 (EP Colour (includes auto selection), USP, Pt-Co)	1371955
PFXi-195/6 (Klett Color (blue filter KS-42), Pt-Co, Hess-Ives)	1371956
PFXi-195/7 (ICUMSA (420, 560, 710 nm), Honey Colour, Series 52)	1371957
PFXi-195/8 (ADMI (spectral & tristimulus filter methods), Pt-Co)	1371958
PFXi-195/9 (full range of CIE values)	1371959

Colour Scale Upgrades (include required cell)

colour state opgrades (include required con)		
Acid Wash Color (ASTM D848)	132950	
ASTM Color	132840	
Chinese Pharmacopoeia	132850	
Gardner Color	132860	
Klett Color (blue filter KS-42); order with instrument only	132930	
Pt-Co/Hazen/APHA scale	132880	
Rosin, US Naval Stores	132890	
Saybolt Colour	132900	
Yellowness Index	132910	
CIE L*C*h	132920	

Many other Scale Upgrades are available. Please contact sales@tintometer.com for details.

The Tintometer Ltd Tel.: +44 (0)1980 664800

sales@tintometer.com

Tintometer - North America

www.lovibondcolour.com

Tel.: +1 941-758-8671

Tintometer - South East Asia

Tel.: +603 3325 2286

Tintometer - China

Spares & Optional Accessories

Tube, 33 mm diameter for ASTM Color

Dark sample cell for CIE Values, 1 mm path length

Flow through control software (for use with PC)

* For use with holder for standard spectrophotometer cells

EBC Colour, set of 5 filters (4, 9.5, 15, 18, 25)

FAC Color, set of 5 filters (7, 13, 15, 29, 39)

EP Brown/Yellow, set of 2 filters (BY3, BY5)

EP Green/Yellow, set of 2 filters (GY2, GY5)

US Pharmacopoeia, set of 3 filters (G, H, P)

User specified filter (specify scale and value)

Certified Reference Materials, 500 ml

ASTM Color < 0.5, (certified to ISO/IEC 17025)

ASTM Color 1, (certified to ISO/IEC 17025)

ASTM Color 3, (certified to ISO/IEC 17025)

ASTM Color 5, (certified to ISO/IEC 17025)

ASTM Color 7, (certified to ISO/IEC 17025)

Gardner Color 2, (certified to ISO/IEC 17025)

Gardner Color 5, (certified to ISO/IEC 17025)

Gardner Color 8, (certified to ISO/IEC 17025)

Pt-Co 500 mg Pt/l (certified to ISO/IEC 17025)

Saybolt Color -10, (certified to ISO/IEC 17025)

Saybolt Color +12, (certified to ISO/IEC 17025)

Saybolt Color +25, (certified to ISO/IEC 17025)

Service pack for RCMSi

Saybolt Color 0, (certified to ISO/IEC 17025)

Pt-Co 5 mg Pt/l, certified under ISO 9001 quality system

Pt-Co 10 mg Pt/l, certified under ISO 9001 quality system

Pt-Co 15 mg Pt/l, certified under ISO 9001 quality system

Pt-Co 30 mg Pt/l, certified under ISO 9001 quality system

Pt-Co 50 mg Pt/l, certified under ISO 9001 quality system

Pt-Co 100 mg Pt/l, certified under ISO 9001 quality system

Flow through cell, 50 mm path length *

Flow through cell, 10 mm path length *

PFXi-195 IQ/OQ Documentation

EP Red, set of 2 filters (R2, R6)

EP Yellow, set of 2 filters (Y2, Y4)

EP Brown, set of 2 filters (B2, B5)

(nominal values quoted)

(W600/OG/10 mm cell with a 9 mm spacer)

Holder for standard spectrophotometer cells (12.5 mm width)

Plastic, disposable cells, 10 mm path length, pack of 100 *

Plastic, disposable cells, 33 mm path length, pack of 50 *

Plastic, disposable cells, 50 mm path length, pack of 50 *

Conformance Filters (nominal values quoted)

Honey Colour (Pfund Equivalents), set of 5 filters (15, 30, 60, 85, 100) 139370 Klett (blue filter KS-42), set of 5 filters (21, 66, 162, 318, 616)

Single filter, certificated (select scale & nominal value from above)

PFXi-195 Lamp 5V

Acid Wash Tube

Replacement sample chamber

Adapter for 10.65 mm tubes

Adapter for 33 mm tubes

Tube, 10.65 mm diameter

Tel.: +86 10 85251111 ext.330





138180

133830

131907

132800

132810

352290

352300

352220

132820

132770

352091

352101

137690

658010

658020

137680

139400

139700

139410

139420

139430

139440

139450

139710

139460

109970

109980

134290

134000

134010

134020

134030

134200

134210

134220

134140

134150

134160

134170

134180

134190

462803

134040

134050

134060

134070

137380