









## Basics to determine the ideal volume of culture medium for use in

### TPP tissue culture vessels

Volume recommendations must be seen as guideline and do not represent an absolute value. TPP is not responsible for any outcome in this regard. The usual media volume to surface area is 0.2-0.5ml/cm<sup>2</sup>, the upper limit is set by gaseous diffusion through the liquid layer and the optimum will depend upon the oxygen requirement of the cells. Also more confluent cultures require more media and changing more frequently.

The number of cells on a confluent plate, dish, or flask will vary with cell type. For this table, we give no recommendations, but indicate source of information such as: Freshny R. Culture of animal cells Wiley-Liss 1994, Google, etc.

#### 1. Tissue Culture Test Plates

Product No.	Version	No. of Wells	Growth area (cm <sup>2</sup> )	Medium recom. (ml)	Medium max. (ml)
92006 92406		6	9.026	3 - 5	10
92012 92412		12	3.466	1 - 2	5
92024 92424		24	1.864	0.5 - 1	3
92048 92448		48	0.875	0.19 – 0.29	0.2
92096 92696		96	0.322	0.1 – 0.25	0.4
92097 92697		96	0.636	0.1 – 0.25	0.4

#### 2. Tissue Culture Dishes

Product No.	Inner-Ø (mm <sup>2</sup> )	Growth area (cm <sup>2</sup> )	Medium recom. (ml)	Medium max. (ml)
93040	34	9.2	2	4
93060	53	22.1	3	15
93100	87	60.0	10	60
93150	137	147.0	20	155

#### 3. Tissue Culture Flasks

Product No.	Version	Growth area (cm <sup>2</sup> )	Medium recom. (ml)	Medium max. (ml)
90025	VENT	25	3 - 5	15
90026	Filter			
90075	VENT	75	8 - 15	65
90076	Filter			
90150	VENT	150	15 - 30	165
90151	Filter			
90300	VENT	300	30 - 50	410
90301	Filter			

#### 4. Clipmax

Product No.	Version	Growth area (cm <sup>2</sup> )	Medium recom. (ml)	Medium max. (ml)
70010	Filter	10	4 - 6	7



### 5. Tissue Culture Tubes

Product No.	Version	Growth area (cm <sup>2</sup> )	Medium recom. (ml)	Medium max. (ml)
91243	Filter	10	2	18
91106	VENT	20	3	13

### 6. Bioreactor (Suspension cells)

Product No.	Max Volume (ml)	Shaker $\varnothing$ mm / speed rpm	Medium recom. (ml)	Medium max. (ml)
86050	50	50 / 250	35	50
87015	15	50 / 250	10	15
87017	15	50 / 180	10	15
87050	50	50 / 150	35	50
87450	450	50 / 150	250	450
87600	600	50 / 150	400	600