## ODO Cap H2

Specification Sheet (Part/REF # 243505)



The standard ODO Cap H0 / H3 is the default option for most applications. The ODO Cap H2 / H4 has a hygienic convex design to minimize bubble accumulation and a PTFE membrane that increases the chemical resistance.

## **Product Specifications**

Wetted Parts	Stainless Steel 1.4435 PTFE (Polytetrafluoroethylene) EPDM (Ethylene propylene elastomer) See compliance details in Material Specification document
Surface Quality of Steel	Ra < 0.4 μm (N5)
Serial Number	Yes
Response Time	t98% < 60 s at 25 °C
Autoclavable	Yes, max. Temperature 130 °C
Steam Sterilizable	Yes, max. Temperature 140 °C
Operating Temperature Range	-10 to 140 °C

Specifications are subject to change without notice https://www.hamiltoncompany.com/sensors/243505 Spec. Version C