

Automated 24/7 Evaporation



MADE IN
GERMANY

The average operational **lifespan of 10 years** is backed by a **3 year warranty** and makes your purchase a worthwhile investment.

Automated in unattended distillation in **two different modes**


Automatic refill of the evaporation flask and continuous **automatic discharge of condensate**

Automatic drainage of residue: residue in the evaporation flask will be drained off

Compatible with all Hei-VAP rotary evaporators (Distimatic Bench-top) and Hei-VAP Industrial large-scale rotary evaporators (Distimatic Industrial)

The automatic module Distimatic Bench-top can be **used in common fume hoods** due to the space saving design

Applications are completed in an automated and unattended process of **solvent separations of nearly 100 %**



Bumping is eliminated by continuously feeding small volumes

Wide range of possibilities to **adjust settings according to your individual needs**

In case of an emergency **the complete system will power off entirely**