

## UNIGLOVES SELECT® PLUS



SIZE	ITEM-NO
------	---------


XS	2301
----	------

S	2302
---	------

M	2303
---	------

L	2304
---	------

XL	2305
----	------



## PRODUCT DESCRIPTION

### Unigloves SELECT® PLUS

#### The grippy glove

The Select Plus is an extremely **non-slip, food-safe** disposable glove made from high-quality natural latex. The hygienic, protective glove is particularly suitable for use in hospitals, medical practices and industry. The Select Plus has been tested in line with **PPE Regulation (EU) 2016/425** and European standards **EN 374, EN 420 and EN 455**. The extremely durable material is waterproof and protects against viruses, fungi and bacteria. The **micro-roughened surface** gives the glove **outstanding grip** and ensures an **excellent hold** when handling smooth or damp objects. This non-powdered glove is easy and quick to put on and remove thanks to chlorination and its rolled cuff. The flexible fit and ambidextrous design allow for fast, frequent changes.

The Select Plus is available in a **natural colour** and sizes **XS, S, M, L and XL**.

## EIGENSCHAFTEN

<b>MATERIAL</b>	<u>Latex</u>
<b>COLOUR</b>	<u>natural latex</u>
<b>LÄNGE</b>	<u>240 mm</u>
<b>SIZE</b>	<u>XS (5.5-6.0), S (6.5-7.0), M (7.5-8.0), L (8.5-9), XL (9.5-10)</u>
<b>PRODUCT FEATURES</b>	<u>beaded cuff, chlorinated on the outside, colour coding, micro-roughened palm, non-sterile, polymer coating on the inside, powder-free, very grippy</u>
<b>AQL</b>	<u>1.5</u>
<b>NORM GLOVES</b>	<u>Lebensmittelecht, CE, EN 374-1 Type-C, EN 455, T. 1-4</u>
<b>VERPACKUNGSINHALT</b>	<u>box of 100 pieces</u>
<b>VERSANDEINHEIT</b>	<u>10 boxes</u>