





UNIGLOVES FORMAT® WHITE

SIZE	ITEM-NO.
XS	1001
S	1002
M	1003
L	1004
XL	1005



PRODUCT DESCRIPTION

Unigloves FORMAT® WHITE

This medical examination glove made of nitrile is primarily used in hospitals and medical practices as well as beauty and hair salons. The disposable glove has been tested in line with PPE Regulation (EU) 2016/425 and European standards EN 374, EN 420 and EN 455. The durable, elastic material is waterproof and protects against viruses, fungi and bacteria. The latex-free, thiuram-free and protein-free glove is particularly skin-friendly. The micro-roughened fingers and palms ensure excellent grip on various surfaces. This non-powdered glove is easy to put on and remove thanks to its strong, cuffed sleeve. The flexible fit and ambidextrous design makes them easy to put on and remove.

The glove is available in white and sizes XS, S, M, L and XL.

EIGENSCHAFTEN

MATERIAL	<i>Nitrile</i>
COLOUR	<i>white</i>
LÄNGE	<i>240 mm</i>
SIZE	<i>XS (5.5-6.0), S (6.5-7.0), M (7.5-8.0), L (8.5-9), XL (9.5-10)</i>
PRODUCT FEATURES	<i>beaded cuff, chlorinated on the inside & outside, colour coding, micro-roughened palm, non-sterile, powder-free</i>
AQL	<i>1.5</i>
NORM GLOVES	<i>CE, EN 374-1 Type-C, EN 455, T. 1-4</i>
VERPACKUNGSINHALT	<i>box of 100 pieces</i>
VERSANDEINHEIT	<i>10 boxes</i>