

Tip compatibility chart

Pipette brand	Pipette model	IKA Tip 10 µl bag Ident. No. 0020114634	IKA Tip 200 µl bag Ident. No. 0020114635	IKA Tip 1000 µl bag Ident. No. 0020114636
		IKA Tip 10 µl box Ident. No. 0020116857	IKA Tip 200 µl box Ident. No. 0020116858	IKA Tip 1000 µl box Ident. No. 0020116859
IKA	PETTE fix single-channel, fixed volume			
	5 µl	x		
	10 µl	x		
	20 µl		x	
	25 µl		x	
	50 µl		x	
	100 µl		x	
	200 µl		x	
	250 µl			x
	500 µl			x
	1000 µl			x
	PETTE vario single-channel, variable volume			
	0.1 – 2 µl	x		
	0.5 – 10 µl	x		
	2 – 20 µl		x	
	10 – 100 µl		x	
	20 – 200 µl		x	
	100 – 1000 µl			x
	0.5 – 5 ml			
	1 – 10 ml			
	PETTE multi 8-channel, variable volume			
	8 x 0.5 – 10 µl	x		
	8 x 10 – 100 µl		x	
	8 x 20 – 200 µl		x	
	PETTE multi 12-channel, variable volume			
	12 x 0.5 – 10 µl	x		
	12 x 10 – 100 µl		x	
	12 x 20 – 200 µl		x	

Tip compatibility chart

Pipette brand	Pipette model	IKA Tip 10 µl bag Ident. No. 0020114634	IKA Tip 200 µl bag Ident. No. 0020114635	IKA Tip 1000 µl bag Ident. No. 0020114636	
		IKA Tip 10 µl box Ident. No. 0020116857	IKA Tip 200 µl box Ident. No. 0020116858	IKA Tip 1000 µl box Ident. No. 0020116859	
EPPENDORF	Reference				
		0.5 – 10 µl	x		
		10 – 100 µl		x	
		50 – 200 µl		x	
		200 – 1000 µl			x
	Research				
		0.5 – 10 µl	x		
		2 – 20 µl		x	
		10 – 100 µl		x	
		20 – 200 µl		x	
		100 – 1000 µl			x
	8-channel	0.5 – 10 µl		x	
	Reference 2 single-channel, variable volume				
		0.5 – 10 µl	x		
		2 – 20 µl		x	
		10 – 100 µl		x	
		20 – 200 µl		x	
		100 – 1000 µl			x
		0.25 – 2.5 ml			
		0.5 – 5 ml			
		1 – 10 ml			
	8-channel	10 – 100 µl		x	
	Research Plus single-channel, variable volume				
		0.5 – 10 µl	x		
	2 – 20 µl	x			
	10 – 100 µl		x		
	20 – 200 µl		x		
	100 – 1000 µl			x	
	0.5 – 5 ml				
	1 – 10 ml				
8-channel	0.5 – 10 µl	x			
8-channel	10 - 100 µl		x		
12-channel	0.5 – 10 µl	x			
12-channel	10 – 100 µl		x		

Tip compatibility chart

Pipette brand	Pipette model	IKA Tip 10 µl bag Ident. No. 0020114634	IKA Tip 200 µl bag Ident. No. 0020114635	IKA Tip 1000 µl bag Ident. No. 0020114636	
		IKA Tip 10 µl box Ident. No. 0020116857	IKA Tip 200 µl box Ident. No. 0020116858	IKA Tip 1000 µl box Ident. No. 0020116859	
SARTORIUS	Proline				
		0.5 – 10 µl	x		
		10 – 100 µl		x	
		20 – 200 µl		x	
		100 – 1000 µl			x
		1000 – 5000 µl			
	mLine				
		0.5 – 10 µl	x		
		10 – 100 µl		x	
		20 – 200 µl		x	
		100 – 1000 µl			x
		500 – 5000 µl			
	GILSON	Pipetman Classic / G / L - traditional shaft			
		P10	1 – 10 µl	x	
P20		2 – 20 µl	x		
P100		10 – 100 µl		x	
P200		50 – 200 µl		x	
P1000		100 – 1000 µl			x
P5000		1 – 5 ml			
P10ML		1 – 10 ml			
Pipetman Ultra					
U10		1 – 10 µl	Ø*		
U100		10 – 100 µl		x	
U200		20 – 200 µl		x	
U1000		200 – 1000 µl			x
U5000		1 – 5 ml			
U10000	1 – 10 ml				

* Tip ejection is not possible.
It has to be manually removed.

Tip compatibility chart

Pipette brand	Pipette model	IKA Tip 10 µl bag Ident. No. 0020114634	IKA Tip 200 µl bag Ident. No. 0020114635	IKA Tip 1000 µl bag Ident. No. 0020114636
		IKA Tip 10 µl box Ident. No. 0020116857	IKA Tip 200 µl box Ident. No. 0020116858	IKA Tip 1000 µl box Ident. No. 0020116859
GILSON	Pipetman Neo			
	1 – 10 µl	x		
	2 – 20 µl	ø*		
	10 – 100 µl		x	
	20 – 200 µl		x	
	100 – 1000 µl			x
SOCOREX	Acura 825 single channel digital			
	0.5 – 10 µl	x		
	2 – 20 µl			
	10 – 100 µl		x	
	20 – 200 µl		x	
	100 – 1000 µl			
	Acura manual 835			
	0.2 – 2 ml			
	0.5 – 5 ml			
	1 – 10 ml			
	Acura manual XS			
	0.5 – 10 µl	x		
	2 – 20 µl			
	10 – 100 µl		x	
	20 – 200 µl		x	
	100 – 1000 µl			x
	Calibra Digital 822			
	1 – 10 µl	x		
2 – 20 µl				
10 – 100 µl				
20 – 200 µl		x		
100 – 1000 µl			x	
Calibra Digital 832				
0.2 – 2 ml				
1 – 10 ml				

Tip compatibility chart

Pipette brand	Pipette model	IKA Tip 10 µl bag Ident. No. 0020114634	IKA Tip 200 µl bag Ident. No. 0020114635	IKA Tip 1000 µl bag Ident. No. 0020114636	
		IKA Tip 10 µl box Ident. No. 0020116857	IKA Tip 200 µl box Ident. No. 0020116858	IKA Tip 1000 µl box Ident. No. 0020116859	
BRAND	Transferpette S				
	0.5 - 10 µl	x			
	2 - 20 µl		x		
	10 - 100 µl		x		
	20 - 200 µl		x		
	100 - 1000 µl			x	
	0.5 - 5 ml				
	1 - 10 ml				
		Transferpette electronic			
	0.5 - 10 µl	x			
	2 - 20 µl		x		
	10 - 100 µl		x		
	20 - 200 µl		x		
	100 - 1000 µl			x	
200 - 5000 µl					
THERMO FINNPIPETTE	Finnpipette digital pipettes and F1 / F2 / F3				
	0.2 - 2 µl	x			
	0.5 - 10 µl	x			
	10 - 100 µl		x		
	20 - 200 µl		x		
	200 - 1000 µl			x	
	1 - 5 ml				
2 - 10 ml					

Please note that the above compatibility table is for guidance only. IKA is not liable for technical modifications or changes to the listed pipettes. The specifications (= precision and accuracy) have not been tested.