

# FT900 Immersion Cooler

## with immersion probe for rapid cooling of liquids

JULABO immersion coolers are primarily employed for countercooling open bath circulators to ambient temperature or below as well as for rapidly cooling liquids. JULABO immersion coolers have a small footprint and represent an economic solution preventing tap water usage for cooling. Ease of operation make this model series the i deal instruments for a variety of cooling applications.

# Your advantages

- · Rapid cooling of liquids down to low temperatures
- · Counter-cooling of liquids in combination with heating circulators
- · Environmentally friendly by saving precious tap water
- · User-friendly operation and handling
- · Compact design, small footprint
- · Dry ice substitution



### **Technical Data**

Order No.	9650890
Model series	FT Series
Category	Immersion Coolers
Working temperature range (°C)	-90 <b>+</b> 30
Cooling capacity (Medium Ethanol)	°C 20 10 -20 -40 -80 kW 0.3 0.27 0.24 0.2 0.07
Refrigerant	R23, R404A
Ambient temperature	535 °C
Dimensions W x L x H (cm)	38 x 55 x 60
Weight (kg)	50
Cooling of compressor	Air
Immersion probe	flexible probe, 65 x 1.5 (L x Ø) cm
Connection tube (L) cm	160
Available voltage versions	230 V / 50-60 Hz 115 V / 60 Hz

JULABO GmbH Eisenbahnstr. 45

77960 Seelbach / Germany

Tel. +49 (0) 7823 51-0 info.de@julabo.com

## **JULABO Services**

Product finder, Accessory search, 1 PLUS Warranty, Catalog download, Callback service, Operating manuals, Safety data sheets, Software and more

www.julabo.com