

## F250 Recirculating Cooler

### for simple cooling tasks

JULABO F models require very little space and have very low procurement costs. Recirculating coolers of the F Series are a great way to replace costly tap water and are ideal for basic cooling tasks.

### Your advantages

- Environmentally-friendly operation with low energy consumption
- Compact design
- Splash-proof membrane keypad with LED temperature display
- Straightforward filling and draining
- Filling level indicator
- May be used with water, water/glycol, JULABO Thermal G



### Technical Data

Order No.	9620025					
Model series	F Series					
Category	Recirculating Coolers					
Working temperature range (°C)	-10 ... +40					
Temperature stability (°C)	±0.5					
Temperature Control	PID temperature control					
Setting / display resolution	0.1 °C					
Temperature Display	LED					
Cooling capacity (Medium Ethanol)	°C	20	10	5	0	-5
	kW	0.25	0.22	0.21	0.18	0.09
Pump capacity flow rate (l/min)	15					
Pump capacity flow pressure (bar)	0.35					
Pump connections	M10x1					
Barbed fittings diameter (inner dia. / mm)	8 / 10					
Filling volume liters	1.7 ... 2.6					
Refrigerant	R134a					
Ambient temperature	5...40 °C					
Dimensions W x L x H (cm)	24 x 40 x 52					
Weight (kg)	27					
Included with each unit	2 each barbed fittings for tubing 8 and 10 mm inner dia. (pump connections M10x1 male)					
Cooling of compressor	Air					
Available voltage versions	230 V / 50 Hz					
	230 V / 60 Hz					

115 V / 60 Hz  
200 V / 50-60 Hz  
100 V / 50-60 Hz

## Characteristics

### Display



#### Easy to read

Large LED temperature display for actual value and setpoint (resolution 0.1 °C)

### Operation



#### Simple and fast

Convenient 3-key setpoint adjustment (F models)

### Temperature Control



#### Precise

PID Temperature control with set control parameters, temperature stability  $\pm 0.02 \dots \pm 0.2$  °C

#### JULABO Contact

JULABO GmbH  
Eisenbahnstr. 45  
77960 Seelbach / Germany

Tel. +49 (0) 7823 51-0  
info.de@julabo.com

#### JULABO Services

Product finder, Accessory search, 1 PLUS Warranty, Catalog download, Callback service, Operating manuals, Safety data sheets, Software and more

[www.julabo.com](http://www.julabo.com)