



810-SG Series **COMPACT GLOVE BOX** with transfer chamber AND flat side access door

Completely sealed glove box for basic isolation and/or containment.





The 810-SG Series by PLAS—LABS are constructed of an optically clear material, developed for isolation and / or containment. All units feature a high visibility of the interior, cleanliness, safety, complete containment, portability and a draft-free atmosphere. The installed purging valves are ideal for lowering oxygen and humidity levels.

Other features include:

- Transfer chamber with purging valves and vacuum gauge.
- Flat side access door for easily introducing larger equipment.
- Four (4) key cock valves for purging: Two on the main body, two on the transfer chamber.

- Hospital grade multiple electrical outlet strip.
- Pressure relief valve with small HEPA filter.
- Standard 24 month warranty, (not on gloves or consumables.)
- White ambidextrous gloves.

	810-SG10	810-SG20	810-SG30	810-SG40
Inside Dimensions	28"w x 23"d x 29"h	35"w x 29"d x 30"h	48"w x 29"d x 32"h	60"w x 29"d x 32"h
	710mm x 580mm x	890mm x 740mm x	1220mm x 740mm x	1524mm x 740mm x
	740mm	760mm	810mm	810mm
Outside Dimensions	43"w x 24"d x 31"h	49"w x 30"d x 31"h	63"w x 31"d x 35"h	76"w x 32"d x 35"h
(Includes transfer chamber 12"	1100mm x 610mm x	1250mm x 760mm x	1600mm x 790mm x	1930mm x 812mm x
long)	790mm	790mm	860mm	890mm
Access Door Opening	15.5" x 22"h	21.5" x 22"h	21.5" x 24"h	21.5" x 24"h
	400mm x 560mm	546mm x 560mm	546mm x 610mm	546mm x 610mm