

BagMixer® 400 P


400 mL lab blender

Ref. 021 000 (100-127 V~ / 220-240 V~)

Ref. 021 230 (220-240 V~)



Standards:



Manufactured under:



BagMixer® 400 P efficiently blends all types of samples in all types of bags, without any risk of contamination. The homogenized sample is representative of the global contamination and ready for analysis.

TECHNICAL SPECIFICATIONS

- All stainless steel 304 L
- Plain door
- Fully opening door (270°)
- Q-Tight® closing: tight closing for no more leaks from the top of the bag
- Also available in CoolMix® version

- Useful volume: 50 - 400 mL
- Fixed speed (8 strokes/s)
- Variable blending time (30 - 210 s or ∞)

- Voltage and frequency
Ref. 021 000: 100-127 V~ / 220-240 V~ 50-60 Hz
Ref. 021 230: 220-240 V~ 50-60 Hz
- Fuses
Ref. 021 000: T 3,15 A L 250 V / T 1,25 A L 250 V 5x20 mm
Ref. 021 230: T 1,25 A L 250 V 5x20 mm
- Max power: 135 W

- Protection index : 21
- Dimensions (w x d x h): 26.8 x 44.1 x 28.3 cm, weight: 16.8 kg
- Box dimensions (w x d x h): 49 x 34.5 x 37 cm, gross weight: 19.5 kg
- 3-year warranty
- Shock absorbers lifetime warranty
- Spare parts availability: 10 years

- In compliance with: ISO 7218, ISO 6887 and FDA-BAM (Bacteriological Analytical Manual)
- Manufactured under CE, RoHS, WEEE
- Designed and made in France

DELIVERED WITH

- 1 power cable
- 1 user's manual
- 1 certificate of conformity / warranty card
- 1 free pack of blender bags

OPTIONS / ACCESSORIES

- Ref 026 000 Drip tray: Stainless steel security drip tray for BagMixer® 400

PARIS

Phone: +33 1 34 62 62 61
info@interscience.com

FRANKFURT

Phone: +49 611 7238 7770
sales.germany@interscience.com

BOSTON

Phone: +1 781 937 0007
sales.usa@intersciencelab.com

SHANGHAI

电话: +86 178 2123 6642
sales.china@interscience.cn

SINGAPORE

Phone: +65 6977 7232
sales.asia@interscience.com

TOKYO

Phone: +81 3 6712 9715
sales.japan@interscience.com