

800 Series MICROLITER™ Syringes

5 µL - 250 µL

- For use with liquids
- Removable needles (RN) or Cemented needles (N)
- Reinforced plungers
- Barrel/plunger assemblies are replaceable
- Plunger stop prevents plunger blowout



N (Cemented Needle)



RN (Removable Needle)




Point 2

Hamilton 800 Series Microliter™ Syringes

Cat.No	Model/Needle Term.	Volume	Needle Length	Needle Gauge	Point Style
HA-84850	85 N	5 µL	51 mm	26s	2
HA-84852	801 N	10 µL	51 mm	26s	2
HA-84854	802 N	25 µL	51 mm	22s	2
HA-84856	805 N	50 µL	51 mm	22s	2
HA-84858	810 N	100 µL	51 mm	22s	2
HA-84860	825 N	250 µL	51 mm	22s	2
HA-84851	85 RN	5 µL	51 mm	26s	2
HA-84853	801 RN	10 µL	51 mm	26s	2
HA-84855	802 RN	25 µL	51 mm	22s	2
HA-84857	805 RN	50 µL	51 mm	22s	2
HA-84859	810 RN	100 µL	51 mm	22s	2
HA-84861	825 RN	250 µL	51 mm	22s	2

Operating Parameters

- Maximum temperature for Cemented Needle is 50°C and Removable Needle 115°C.
- Maximum test pressure is 6 bar



Replacement Parts

Description	Volume Model	5 µL	10 µL	25 µL	50 µL	100 µL	250 µL
		85	801	802	805	810	825
1. Barrel/Plunger Assembly N		HA-32150	HA-32151	HA-32152	HA-32153	HA-32154	HA-32155
Barrel/Plunger Assembly SN ¹⁾		HA-32105	HA-32158	HA-32159	HA-32160	HA-32161	HA-32162
Barrel/Plunger Assembly RN		HA-32164	HA-32165	HA-32166	HA-32167	HA-32168	HA-32169
B/P Assembly Syringe ²⁾		HA-32134	HA-32129	HA-32117	HA-32120	HA-32123	HA-32126
2. Holder Assembly		HA-32135	HA-32135	HA-32135	HA-32135	HA-32135	HA-32135
3. RN Nut		HA-30902	HA-30902	HA-30902	HA-30902	HA-30902	HA-30902
4. RN Needle See specific page							

¹⁾ Specify needle gauge, length, point style, and whether electro-trapped

²⁾ Needle not included