



Product description

The new LLG-proMLP single channel microliter pipettes combine high accuracy and reliability with ergonomic and extremely low force pipetting relieving stress from your hand and thumb. This high precision instrument has been built to maintain accuracy for many years making it the ideal choice for all scientific work requiring exact and reproducible pipetting. The LLG-proMLP pipettes are equipped with an adjustable eject button allowing right and left handed users to set the button* according to their preferences. The built in volume locking mechanism and competitive pricing complete these high quality precision pipettes. Available volume ranges from 0.2 µl to 10 ml.

*not available for 5 ml and 10 ml pipettes

Characteristics

- High quality Made in France
- High accuracy, reliability and reproducibility
- Extremely light weight while the robustness still remains
- Perfect ergonomics for highest comfort
- Colour-coded push-button for quick operation
- Extremely low pipetting forces
- Built-in volume locking for error-free results
- Tip ejector button, adjustable for right- and left-handed users
- Completely autoclavable at 121 °C
- Including individual calibration certificate (according to EN ISO 8655)

Type	Capacity µl	Increments µl	Accuracy max. vol. (≤ VK%)	Präcision max. vol. (≤ CK%)	Colour code	PK	Cat. No.
LLG-proMLP	0,2-2	0.002	1.65	0.75	dark orange	1	6.270 420
LLG-proMLP	0,5-10	0.02	1.10	0.44	red	1	6.270 421
LLG-proMLP	2-20	0.02	1.00	0.33	light yellow	1	6.270 422
LLG-proMLP	10-100	0.2	0.80	0.17	light orange	1	6.270 423
LLG-proMLP	20-200	0.2	0.80	0.17	dark yellow	1	6.270 424
LLG-proMLP	100-1000	2	0.80	0.17	dark blue	1	6.270 425
LLG-proMLP	500-5000	2	0.66	0.18	purple	1	6.270 426
LLG-proMLP	100-10000	20	0.60	0.18	light blue	1	6.270 427