





The new RV 3 Rotary Evaporator is the ideal entry-level model of the IKA rotary evaporator portfolio. It finds a multitude of uses in the chemical, pharmaceutical, and biotechnology industries, in research and development, in manufacturing and quality assurance, in laboratories, and in plant construction. Thanks to specially designed glass guides, the condenser makes extremely efficient use of the 1500 cm<sup>2</sup> cooling surface.

- New: 4 I heating bath for up to 100°C
- water heating bath with digital temperature display and carrying handles
- mechanical lift end-point safety stop



- locking mechanism: red indicator shows unlocked position of the vapor tube
- manual lift for precise positioning of the glassware
- adjustable immersion angle
- single-handed manual lift handling, suitable for left and tight-handed operators
- stepless speed setting with dial control and speed display
- speed range: 5 300 min-1
- low device voltage (24V) ensures user safety
- flask clamping mechanism with integrated push-off function for easy exchange of evaporation flasks
- high-efficiency condenser with 1500 cm<sup>2</sup> cooling surface low space requirements
- compatible with the entire range of IKA RV 10 glassware



## **Technical Data**

Type of cooling	vertical	
Cooling surface [cm2]	1500	
Motor principle	DC	
Speed range [rpm]	20 - 300	
Speed tolerance set rotation speed < 100rpm [±rpm]	1	
Speed tolerance set rotation speed > 100rpm [±%]	1	
Lift	manual	
Stroke [mm]	150	
Heating temperature range [°C]	Room temp 90	
Dimensions (W x H x D) [mm]	440 x 530 x 330	
Weight [kg]	13.5	
Permissible ambient temperature [°C]	5 - 40	
Permissible relative humidity [%]	80	
Protection class according to DIN EN 60529	IP 20	
Voltage [V]	100 - 240	
Frequency [Hz]	50/60	
Power input [W]	1400	