



## Product data sheet

**LLG-Disposable Shoe covers**  
**Art.-Nr.: 6.282 795**

Material: CPE  
Weight: 3.3g+/- 2g/pcs  
Color: blue  
Size: 15x41cm

CE MARK  
ISO 13485: 2003



### Safety instructions:

1. Identification of Hazards

No hazards.

2. First-aid measures

In case of rashes or any other discomfort occur, discontinue use immediately and wash hand with plenty of water. Consult physician.

3. Fire-fighting measures

Extinguishing media : Water, foam or dry powder

In case of fire, wear respiratory protection equipment.

4. Handling and storage

Store in cool and dry area.

Protect from direct sun light and UV light.

Avoid prolonged exposure to atmospheric air.

5. Toxicology information

Skin irritation: No irritation (rabbit)

Dermal contact sensitization No sensitization (guinea pig)