

JumboMix® 3500 WarmMix® CC


3500 mL lab blender

Ref. 031 110 (100-127 V~ / 220-240 V~)

Ref. 031 230 (220-240 V~)



Standards:



Manufactured under:



JumboMix® 3500 WarmMix® CC enables the quick homogenizing of large samples. Avoiding all risks of cross-contamination, it is an easy-to-use and powerful lab blender. Its regulated door (from 20°C to 50°C) provides perfect conditions for the Trichinella testing in pork.

TECHNICAL SPECIFICATIONS

- All stainless steel 304 L
- Multifunction digital screen
- Regulated heated door from 20° C to 50° C (68° F to 122° F)
- Fully opening door (270°)
- Q-Tight® closing: tight closing for no more leaks from the top of the bag
- Click & Clear® removable paddles (patented)
- Side-by-side paddle stop
- Security drip tray

- Adjustable blending power (patented)
- Optimal blending volume: 200 - 3750 mL
- Variable speed (1.5 - 3 - 4.5 - 6 strokes/s)
- Variable blending time (1 s - 1 h or ∞ / count-down)

- Voltage and frequency
Ref. 031 110 : 100-127 V~ / 220-240 V~ 50-60 Hz
Ref. 031 230 : 220-240 V~ 50-60 Hz
- Fuses
Ref. 031 110 : T4 A L 250 V / T2 A L 250 V 5x20 mm
Ref. 031 230 : T2A L 250 V 5x20 mm
- Max power: 400 W

- Protection index : 21
- Dimensions (w x d x h): 61.6 x 52.3 x 46.8 cm, weight: 57.8 kg
- Box (w x d x h) : 61 x 57 x 57 cm, weight: 62.5 kg
- 3-year warranty
- Shock absorbers lifetime warranty
- Availability of spare parts : 10 years

interscience

- In compliance with: ISO 7218, ISO 6887-1 and FDA BAM (Bacteriological Analytical Manual)
- Manufactured under CE, RoHS, WEEE - Products made for INTERSCIENCE by Interlab, an ISO 9001 certified company
- Designed and made in France

DELIVERED WITH

- 1 security drip tray
- 1 power cable
- 1 user's manual
- 1 quick user guide
- 1 free pack of blender bags

OPTIONS / ACCESSORIES

- Ref 037 000 JumboStand: Ergonomic table for JumboMix® 3500

PARIS

Phone: +33 1 34 62 62 61
info@interscience.com

FRANKFURT

Phone: +49 611 7238 7770
sales.germany@interscience.com

BOSTON

Phone: +1 781 937 0007
sales.usa@intersciencelab.com

SHANGHAI

电话: +86 178 2123 6642
sales.china@interscience.cn

SINGAPORE

Phone: +65 6977 7232
sales.asia@interscience.com

TOKYO

Phone: +81 3 6712 9715
sales.japan@interscience.com