



LABORHITS

LLG-Safety Shields Series







basic

Perfect for visitors as well as for spectacle wearers



- CE approved according to EN166 and EN170 $\,$
- Perfect for visitors as well as over goggles for wearers of glasses
- Integrated side lens with venting, 180° visual field
- Frame made of tested hypo-allergenic material
- 100 % UV protection (up to 380 nm) Weight: 43 g

Colour	Lens	PK	Cat. No.
clear	clear, uncoated, 2-1.2 U 1 F CE	10	9.006 002

classic

ergonomic fit & attractive price



- CE approved according to EN166 and EN170 $\,$
- Integrated side-lenses, 180° visual field Frame made of tested hypo-allergenic material
- Side arms adjustable in length and inclination
- 100 % UV protection (up to 380 nm)
- Weight: 34 g only

Colour	Lens	PK	Cat. No.
blue	clear, anti-scratch lens, 2C-1.2 U 1 F CE	1	9.006 001

shields Series



classic light

perfect fit





- CE approved according to EN166 and EN170
- Integrated side-lenses, 180° visual field
- Frame made of tested hypo-allergenic material
- Side arms adjustable in length
- perfect fit
- 100% UV protection (up to 380 nm)
- Weight: 29 g only

Colour	Lens	PK	Cat. No.
light blue	clear, AF, AS, 2C-1.2 U 1 FT CE	1	9.006 011

evolution and evolution+

also available tinted



- CE approved according to EN166, EN170 and EN 172 for tinted lens
- Panoramic wrapping lens, 180° visual field
- Frame made of tested hypo-allergenic material
- Side arms adjustable in length
- Nose pad and arm tips made of soft, mouldable material offer optimal comfort and individual customisation
- Rubber eyebrow protection for better shock absorption
- 100 % UV protection (up to 400 nm)
- Weight: 30 g only

_				
Тур	Colour	Lens	PK	Cat. No.
Evolution	black/green	clear, Anti-Scratch +, Anti-Fog +	1	9.006 010
Evolution+	black/green	tinted, Anti-Scratch +, Anti-Fog +	1	9.006 000





LABORHITS


blue comfort & fashion length soft ends adjustable 180°

- CE approved according to EN166 and EN170 $\,$
- Panoramic wrapping lens, 180° visual fiel
- Frame made of tested hypo-allergenic material
- Side arms adjustable in length and inclination
- Soft temple ends made with outsert moulding
- 100 % UV protection (up to 380 nm)
- Weight: 29 g only

Colour	Lens	PK	Cat. No.
blue/silver	Clear, anti-scratch lens, 2C-1.2 U 1 F CE	1	9.006 012

lady

slim fit



- CE approved according to EN166 and EN170
- Panoramic wrapping lens, 180° visual field
 Frame made of tested hypo-allergenic material
- Temples with rubber tips, adjustable in length and inclination
- Nosepad and temple tips made of soft, mouldable material offer optimal comfort and individual customisation
- 100 % UV protection (up to 400 nm)
- Weight 25 g only

Colour	Lens	PK	Cat. No.
white/green	clear, AF+, AS+, 2C-1.2 U 1 FT KN CE	1	9.006 007
white/pink	clear, AF+, AS+, 2C-1.2 U 1 FT KN CE	1	9.006 008
black/green	clear, AF, AS, 2C-1.2 U 1 FT CE	1	9.006 009