

LAB & INDUSTRY EQUIPMENT

BUCH  HOLM



Buch & Holm A/S

Marielundvej 39 • 2730 Herlev • Tlf. +45 44 54 00 00 • e-mail b-h@buch-holm.dk • www.buch-holm.dk
Box 10058 • S-200 43 Malmö • Tel. +46 (0)40 30 50 30 • e-mail b-h@buch-holm.se • www.buch-holm.se

DuraCal pH buffer 4.01/7.00/10.01, 500 ml each

Specification Sheet (Part/REF # 238924)

Patented buffer solutions provide pH stability of up to five years. High buffering capacity provides rapid and stable calibration.

Product Specifications

| | |
|---|-------------------------------|
| pH Accuracy | ± 0.01 / ± 0.02 pH |
| Package | 500 mL, each |
| Stability Months (after date of manufacturing) | ± 0.01 pH: 24 / ± 0.02 pH: 60 |
| Certificate | DAkkS |
| pH Value | 4.01/7.00/10.01 |

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/process-analytics/sensors/238924>