

Thermostatically controlled incubators - TC series

with standard / glass door



The TC series of thermostatically controlled cabinets is used for continuous temperature control over a range of 2 °C to 40 °C. This makes them ideal for a wide range of different applications in industrial and research laboratories.

In particular they are ideal for the temperature-controlled storage of samples or BOD determination in effluent analysis work.

The temperature can be set in steps of 0.1 °C and an LED display shows both the set temperature and the current temperature in the cabinets. Devices such as magnetic agitators, which require a power supply, can be connected to sockets incorporated in the interior of the cabinet. The integral temperature control unit meets the requirements of the EMC directive issued as IEC 61326: "Electrical devices for measurement, monitoring and for use in laboratories".

Improved, robust, insulated housing and highly efficient components provide maximum energy efficiency.

There are 4 models available with standard doors from 135 to 445 litres net capacity, and 2 models with glass doors with 140 and 255 litres net capacity.

Highlights

- Temperature range 2 °C to 40 °C, continuously adjustable in steps of 0.1 °C
- Low power consumption
- Illuminated LED display of preset and current temperatures
- Ideal for BOD determination at 20 °C
- Power sockets inside the incubator
- 6 models in 4 sizes
- Standard door or glass door



Temperature control unit

The temperature control unit fulfills the EMC requirements according IEC 61326 : Electrical equipment for measurement, control and laboratory use.

Technical Data

Design	Fully insulated cabinet with universal temperature control unit
Lock	existing
Models with glass door	Insulating glass door in an ABS frame. Ceiling lighting, separately switchable
Operation	Splash-protected keypad, 2 buttons with tactile feedback
Control range	+ 2 °C to + 40 °C, steps of 0.1 °C
Climate class	+ 10 °C to + 32 °C,
Temperature tolerance	± 1 °C, specified for a stirred 500 ml water sample. For BOD (T=20 °C ±0,5 °C)
Display	Backlit LED display, Resolution 0.1 °C
Fan	Axial, output 320 m ³ /h
Cooling/Heating	Integrated powerful cooling and heating
Power supply	220 - 240 V / 50 Hz
Sockets	CEE 7/5, type E with hinged lid, 230 V / 16 A, 2p + E, IP 44
Coolant	R134a
Approval	CE

Models with standard door

TC 135 S

3 metal racks + 1 bottom grid + 4 sockets

Consumption: approx. 1.35 kWh / 24 h*

I. D. (approx.): 513 W x 441 D x 702 H mm

Net capacity: approx. 135 l

O. D. (approx.):

600 W x 600 D x 850 H mm with work top

600 W x 600 D x 819 H mm without work top

Suitable for built under applications

Weight: approx. 39.0 kg

Order code: 438200



Models with glass door

TC 140 G

3 metal racks + 1 bottom grid + 4 sockets

Consumption: approx. 1.77 kWh / 24 h**

I. D. (approx.): 513 W x 441 D x 702 H mm

Net capacity: approx. 140 l

O. D. (approx.):

600 W x 600 D x 850 H x mm with work top

600 W x 600 D x 819 H mm without work top

Suitable for built under applications

Weight: approx. 48.0 kg

Order code: 438210



TC 175 S

3 metal racks + 1 bottom grid + 5 sockets

Consumption: approx. 1.23 kWh / 24 h*

I. D. (approx.): 470 W x 440 D x 1062 H mm

Net capacity: approx. 175 l

O. D. (approx.): 600 W x 610 D x 1250 H x mm

Weight: approx. 51.0 kg

Order code: 438220



TC 256 G

4 metal racks + 1 bottom grid + 7 sockets

Consumption: approx. 1.56 kWh / 24 h**

I. D. (approx.): 470 W x 440 D x 1452 H mm

Net capacity: approx. 255 l

O. D. (approx.): 600 W x 610 D x 1640 H x mm

Weight: approx. 77.0 kg

Order code: 438235



TC 255 S

4 metal racks + 1 bottom grid + 7 sockets

Consumption: approx. 1.54 kWh / 24 h*

I. D. (approx.): 470 W x 440 D x 1452 H mm

Net capacity: approx. 255 l

O. D. (approx.): 600 W x 610 D x 1640 H x mm

Weight: approx. 61.0 kg

Order code: 438230



** Ambient temperature 25 °C

Target temperature 20 °C

with interior lighting switched on (15 W)

Variations possible

TC 445 S

4 metal racks + 1 bottom grid + 9 sockets

Consumption: approx. 1.42 kWh / 24 h*

I. D. (approx.): 600 W x 560 D x 1452 H mm

Net capacity: approx. 445 l

O. D. (approx.): 750 W x 730 D x 1640 H x mm

Weight: approx. 78.5 kg

Order code: 438240



* Ambient temperature 25 °C

Target temperature 20 °C

Variations possible

Space for BOD BD 600 systems in thermostatically controlled incubators

Model	systems, standard ¹⁾	systems, comfort ²⁾
TC 135 S / TC 140 G	3	2
TC 175 S	5	2
TC 255 S / TC 256 G	7	3
TC 445 S	12	9

¹⁾ Change of bottles **by** removing racks.

²⁾ Change of bottles **without** removing racks.

Applications

- BOD-Measurement
- Microbiological Research
- Food Industry
- Dairies
- Laboratories
- Research Centres
- Universities