

DuraCal pH buffer 4.01/7.00/10.01, 500 ml each

Specification Sheet (Part/REF # 238924)

Patented buffer solutions provide pH stability of up to five years. High buffering capacity provides rapid and stable calibration.

Product Specifications

pH Accuracy	± 0.01 / ± 0.02 pH
Package	500 mL, each
Stability Months (after date of manufacturing)	± 0.01 pH: 24 / ± 0.02 pH: 60
Certificate	DAkkS
pH Value	4.01/7.00/10.01

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/process-analytics/sensors/238924>